Whatman™ 25mm GD/XP Syringe Filters

Product Information sheet

Warning

For research use only.

Not recommended or intended for diagnosis of disease in humans or animals.

Do not use internally or externally in humans or animals.

The Whatman 25mm GD/XP Syringe Filters are designed to prepare and enhance the inorganic ion analysis of hard-to-filter samples greater than 10 mL. They are available in a wide variety of filter choices with a polypropylene housing and utilize advanced construction methods and design features.

Each unit contains a prefiltration stack of polypropylene prefilters. This Prefiltration Technology allows you to filter difficult samples with reduced hand force.

Disposable filtration devices provide labor saving efficiency while ensuring consistent filter performance when compared to hand assembled reusable filter housings.

This document provides general information on the products listed below. The specifications in the Technical Data section are intended to provide a basis for establishing functional use, as well as for setting quality assurance test performance levels.

Whatman 25mm GD/XP Syringe Filters

Whatman 25mm GD/XP Syringe Filters are designed to prepare and enhance the inorganic ion analysis of hard-to-filter samples greater than 10 mL.

Filter Media	Typical Application
dp PP	High flow and loading capacity aqueous and organic samples
NYL	Aqueous and/or organic samples; hydrophilic
PP	Aqueous and organic samples
PES	Aqueous based samples; low protein binding
PTFE	Organic based samples. Hydrophobic membrane
PVDF	Aqueous and/or organic based samples; low protein binding membrane

Polypropylene Depth Prefilter Considerations:

Compatible with aqueous and/or organic samples, the depth prefiltre stack of polypropylene media, gives extremely high loading capacity, fast flow rates and larger sample volume throughout.

Nylon Membrane (NYL) Considerations:

Nylon membrane is hydrophilic and is a good choice for aqueous and/or aqueous organic samples. The membrane offers good chemical resistance to most common HPLC solvents, however it has limited resistance to acids, bases, halogenated hydrocarbons, aldehydes and strong oxidizing agents. The most common application is HPLC sample filtration.

Polyethersulfone (PES) Membrane Considerations:

This carefully selected Polyethersulfone membrane meets the special filtration requirements of durability, high temperature resistance, good chemical compatibility, and low protein adsorption. It is particularly suitable for filtration of serum, plasma and tissue culture solutions as well as other protein containing solutions where minimal adsorptive protein loss is desired.

Polypropylene Depth Filter (dpPP) Considerations:

This highly retentive polypropylene medium gives excellent flow rates and is compatible with both organic and aqueous samples. Together with polypropylene prefiltration stack it provides a pure, single material construction filter device, with a broad range of solvent/chemical resistance.

Polypropylene (PP) Membrane Considerations:

This specially selected Polypropylene membrane provides high flow with both solvent and aqueous compatibility. This medium provides a pure, single-material-construction filter device with a broad range of solvent/chemical resistance.

Polytetrafluoroethylene (PTFE) Membrane Considerations:


Polytetrafluoroethylene membrane is hydrophobic and will not allow water to pass without high pressures. Aqueous solutions may be filtered if the membrane is initially "wetted" with alcohol or another appropriate solvent. Polytetrafluoroethylene membrane will stop aqueous aerosols in gas streams.

Polyvinylidene Fluoride (PVDF) Membrane Considerations:

Polyvinylidene Fluoride membrane is a suitable choice for most HPLC Sample Prep applications. The membrane is slightly hydrophobic with low water breakthrough values. It offers good chemical resistance to common HPLC solvents, has low protein binding and negligible extractables.

Technical Data

Whatman 25mm GD/XP Syringe Filters

Dimensions:

21.7 mm x 29.7 mm

Weight: Approximately 3 grams

4.6 cm²

Maximum Pressure: 75 psi

Polypropylene

Volume "Hold Up": Full housing 1.4 mL

with air purge 250 µL (approx)

Flow Direction: Flow should enter from the inlet

Inlet - Female Luer Lock (FLL)

Outlet - Male Slip Luer (MSL)

Autoclaving: Autoclave at 121°C (131°C max) at 15 psi for 20 minutes
(Not recommended for Nylon)

Filter Media: See chart below

Integrity Test Data

Description	Pore Size Microns	Minimum Bubble Point psi
dp Polypropylene (dp PP)	0.45	2.4
Nylon (NYL)	0.45	29
Polyethersulfone (PES)	0.45	22
Polypropylene* (PP)	0.45	22
Polytetrafluoroethylene* (PTFE)	0.45	8.5
Polyvinylidene Fluoride (PVDF)	0.45	28

*Bubble point determined with IPA. All others determined with water.

Chemical Compatibility of Membrane

Solvent	DpPP	NYL	PES	PP	PTFE**	PVDF
Acetic Acid, 5%	R	R	R	R	R	R
Acetic Acid, Glacial	R	LR	R	R	R	R
Acetone	R	R	NR	R	R	NR
Acetonitrile	R	R	NR	R	R	R
Ammonia, 6N	R	R	R	R	R	LR
Amyl Acetate	R	R	LR	R	R	LR
Amyl Alcohol	R	R	NR	R	R	R
Benzene*	NR	LR	R	NR	R	R
Benzyl Alcohol*	R	LR	NR	R	R	R
Boric Acid	R	LR		R	R	R
Butyl Alcohol	R	R	R	R	R	R
Butyl Chloride*	NR	NR		NR	R	R
Carbon Tetrachloride*	NR	LR	NR	NR	R	R
Chloroform*	LR	NR	NR	LR	R	R
Chlorobenzene*		NR	NR	LR	R	R
Citric Acid		LR	R	R	R	R
Cresol*	NR	NR	NR	NR	R	NR
Cyclohexanone	R	NR	NR	R	R	R
Cyclohexane*	NR	NR	NR	NR	R	R
Diethyl Acetamide	R	R		R	R	NR
Dimethyl Formamide	R	R	NR	R	R	NR
Dioxane	R	R	LR	R	R	LR
DMSO	R	R	NR	R	R	LR
Ethanol	R	R	R	R	R	R
Ethers*	NR	R	R	NR	R	LR
Ethyl Acetate	R	R	NR	R	R	NR
Ethylene Glycol	R	R	R	R	R	R
Formaldehyde*	LR	R	R	LR	R	R
Freon TF*	NR	NR	R	NR	R	R
Formic Acid	R	NR	R	R	R	R
Hexane	R	R	R	R	R	R
Hydrochloric Acid, Conc*	LR	NR	R	LR	R	R
Hydrofluoric Acid*	LR	NR		LR	R	R
Isobutyl Alcohol	R	R		R	R	R
Isopropyl Alcohol	R	R		R	R	R
Methanol	R	R	R	R	R	R
Methyl Ethyl Ketone	R	R	NR	R	R	NR
Methylene Chloride*	LR	NR	NR	LR	R	R
Nitric Acid, Conc*	NR	NR	NR	NR	R	R
Nitric Acid, 6N*	LR	NR	LR	LR	R	R
Nitrobenzene*	R	LR	NR	R	R	R
Pentane*	NR	R	R	NR	R	R
Perchloro Ethylene*	NR	LR	NR	NR	R	R
Phenol 0.5%	R	NR	NR	R	R	R
Pyridine	R	LR	NR	R	R	NR
Sodium Hydroxide, 6N	R	LR	R	R	R	NR
Sulfuric Acid, Conc*	NR	NR	NR	NR	R	NR
Tetrahydrofuran*	LR	R	NR	LR	R	R
Toluene*	LR	LR	NR	LR	R	R
Trichloroethane*	LR	LR	NR	LR	R	R

Solvent	DpPP	NYL	PES	PP	PTFE**	PVDF
Trichloroethylene*	LR	NR	NR	LR	R	R
Water	R	R	R	R	R	R
Xylene*	LR	LR	LR	LR	R	R

Note: (R = Resistant; LR = Limited Resistance; NR = Non Resistant; + = Insufficient Data; * = Short term resistance of housing).

** = membrane may need pre-wetting with isopropanol/methanol if filtering a polar liquid

Operating Instructions

Safety: When considering the special factors of your application, consult the Technical Data to determine correctness of use. Do not exceed the pressure, temperature or chemical compatibility recommendations. High pressures can be obtained when using syringes. The smaller the syringe the higher the pressure that can be generated. As a general guide, the following pressures can be obtained by hand with the syringes indicated: 20 mL, 80 psi; 10 mL, 140 psi; 5 mL, 180 psi; 3 mL, 200 psi; 1 mL, 250 psi. Each user should determine the pressure they can generate by hand with a specific size syringe and take appropriate safety precautions not to exceed the recommended rating for the device used. If these limitations are exceeded, bursting of the device may occur resulting in loss of sample or personal injury.

Pre-Filter Media: The Whatman 25mm GD/XP Syringe Filters contain a prefiltration stack of proprietary polypropylene prefilters. The first prefILTER is rated at approximately 10 µm while the second is rated at approximately 5 µm. Whatman Prefiltration Technology allows you to filter even the most difficult samples with less hand force than ever before possible. Compared to an unprotected membrane, the volume of sample filtered can be three to seven times greater.

Efficiency: To maximize filtration throughput, use the largest pore size filter that will provide the required cleanliness. To extend filter life use low flow rates or pressures.

To Use With a Syringe:

- 1) Fill the syringe with the solution to be filtered.
- 2) Secure the filled syringe to the FLL inlet of the syringe filter with a twisting motion.
- 3) With the outlet pointed upward, gradually apply thumb pressure to the syringe plunger to initiate flow.
- 4) Continue thumb pressure until all the air in the device is displaced with liquid.
- 5) Once liquid starts to exit syringe filter outlet, stop applying pressure, point the device downward and away from user.
- 6) Position syringe filter over suitable collection container or other apparatus and apply pressure again to filter sample. Change filters when flow becomes too slow or resistance becomes excessive.

Air Locks: Seriously hamper flow rates. To eliminate, point the outlet of the filter device upward during the initiation of liquid flow.

Bubble Point (BP) Test: Flush the filter device with 1.0 ml or more of the test fluid. After the filter is completely wet and with the outlet pointed upward, apply air under controlled pressure to the inlet until air breaks through the filter and bubbles can be seen exiting the outlet. The pressure at which air passes through the wetted filter is the BP. Refer to table for typical BP values.

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Note: If there are any questions regarding these products, a variation of these products, or any Whatman products, please contact our Technical Services Department toll free at 800.922.0361.

Legal.

GE, imagination at work and GE Monogram are trademarks of General Electric Company.

Whatman is a trademark of GE Healthcare companies.

© 2007–2013 General Electric Company – All rights reserved.

Previously published Sept 2007.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request.

Contact your local GE Healthcare representative for the most current information.

For your local office contact information, visit

www.gelifesciences.com/contact

GE Healthcare UK Limited

Amersham Place

Little Chalfont, Buckinghamshire,
HP7 9NA, UK

<http://www.gelifesciences.com>

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala,
Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg,
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327,
Piscataway, NJ 08855-1327,
USA

GE Healthcare Japan Corporation
Sanken Bldg. 3-25-1, Hyakunincho,
Shinjuku-ku, Tokyo 169-0073,
Japan

imagination at work