

BCS Assay Kit

Code	Description	Size
N962-500RXN	BCS Assay Kit	500 reactions
N962-100RXN	BCS Assay Kit	100 reactions
N962-SAMPLE	BCS Assay Kit	20 reactions
N962-Q-SAMPLE	BCS Assay Kit	20 reactions

General Information

The BCS Assay Kit allows for spectrophotometric determination of protein concentration. The BCS Assay utilizes a biuret-like reaction where Copper (II) binds protein. The remaining free Copper (II) then is reduced to Copper (I), which associates with a dye that changes colors upon binding of copper (I). Protein concentration is then determined by measuring the absorbance of the dye-Copper (I) complex. The signal is inversely proportional to the concentration of peptide bonds as opposed to specific amino acid side chains. Therefore, protein to protein variability is less of a concern. Another advantage of the BCS Assay is that the BCS Assay has a high tolerance for salts, detergents, and buffers that may interfere with conventional protein quantitation kits.

- Reduced variability between protein assays
- High tolerance for salts, detergents and buffers
- Reproducible and consistent results

Storage/Stability

Store room temperature (18 to 26°C).

Product Use Limitations

For research use only. Not for therapeutic or diagnostic use.

Materials Supplied

Component	N962-500RXN	N962-100RXN	N962-SAMPLE
BCS Assay Kit, Solution A	N959-50ML	N959-10ML	N959-2ML
BCS Assay Kit, Solution B	N960-450ML	N960-90ML	N960-18ML

Protocol/Procedure:

Standard Curve/Sample Preparation

- For optimal results, BSA standards should be prepared in the final storage buffer of the unknown protein preparation.
- Prepare BSA standards according to the table below.

Standard	Diluent (μL)	20 mg/mL BSA (μL)	Final Concentration (mg/mL)
1	400	0	0
2	396	4	0.2
3	392	8	0.4
4	388	12	0.6
5	384	16	0.8
6	380	20	1

- Aliquot 50 μL of each standard into a 1.5 mL microcentrifuge tube
Note: To increase accuracy, standard curve and samples should be aliquoted and measured in triplicate.
- Aliquot 50 μL of the unknown sample into 1.5 mL microcentrifuge tubes (in triplicate).

Protein Detection

- Add 100 μL of BCS Assay Kit, Solution A to each sample
- Briefly vortex to mix
- Add 900 μL of BCS Assay Kit, Solution B to each sample
- Briefly vortex to mix
- Measure the absorbance of each sample at 485 nm against water or a diluent reference.
- Construct a standard curve by plotting the absorbance of the BSA standards against their respective concentrations. From this graph, the concentration of the unknown protein preparation can be calculated.




Figure. Representative BSA standard curve. The error bars represent the standard deviation obtained from absorbance of three separate samples for each concentration.

Compound	Concentration
Glycerol	≤ 10% (v/v)
Triton X-100	≤ 1% (v/v)
Tween-20	< 0.2% (v/v)
Nonidet ®P-40 Substitute	≤ 2% (v/v)
EDTA	≤ 0.1 mM
DTT	≤ 0.3 mM

Table. BCS Assay reagent compatibility. For protein preparation storage buffers containing reagents above the concentrations mentioned, the BSA standards should be prepared with storage buffer. Alternatively, the unknown protein preparation can be diluted to minimize the effects of the reagents mentioned in the table.



For Technical Support

Toll Free: 1-800-610-2789 (USA & Canada)

Fax: (440) 349-0235

Email: techinquiry@amresco-inc.com

AMRESKO, LLC

A VWR Company

Corporate Headquarters

28600 Fountain Parkway

Solon, Ohio USA 44139-4300

Tel: 440/349-1199

Fax: 440/349-1182

www.amresco-inc.com

BCS Assay Kit

ZY0612

Rev. 1 12/2015

© Copyright 2010 by AMRESKO, LLC

All Rights Reserved.

AMRESKO® is a registered trademark of AMRESKO, LLC