

Product Information

CF™350 Conjugated Antibodies

Catalog No.	Unit Size	Product Description
20140	0.5 mL	Goat Anti-Mouse IgG (H+L)
20140-1	50 uL	
20141	0.5 mL	Goat Anti-Rabbit IgG (H+L)
20141-1	50 uL	
20142	0.5 mL	Donkey Anti-Goat IgG (H+L), highly cross-adsorbed (min X Chicken, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, and Syrian Hamster)
20142-1	50 uL	
20143	0.5 mL	Goat Anti-Mouse IgG (H+L), highly cross-adsorbed (min X Bovine, Horse, Human, Rabbit, and Swine)
20143-1	50 uL	
20144	0.5 mL	Goat Anti-Rabbit IgG (H+L), highly cross-adsorbed (min X Human, Mouse, and Rat)
20144-1	50 uL	
20145	0.25 mL	Goat Anti-Mouse IgG (H+L) F(ab') ₂ fragment
20146	0.25 mL	Goat Anti-Rabbit IgG (H+L) F(ab') ₂ fragment
20147	0.25 mL	Goat Anti-Rat IgG(H+L), highly cross-adsorbed (min X Bovine, Horse, Human, and Rabbit)
20147-1	0.5 mL	
20148	50 uL	Donkey Anti-Sheep IgG (H+L), highly cross-adsorbed (min X Chicken, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, and Syrian Hamster)
20148-1	0.5 mL	
20149	0.5 mL	Rabbit Anti-Mouse IgG (H+L), highly cross-adsorbed (Min X Human)
20149-1	50 uL	
20198	0.5 mL	Goat Anti-Guinea Pig IgG (H+L)
20198-1	50 uL	
20245	0.25 mL	Goat anti-mouse IgG1 (γ1)
20255	0.25 mL	Goat anti-mouse IgG2a (γ2a)
20265	0.25 mL	Goat anti-mouse IgG2b (γ2b)
20275	0.5 mL	Donkey Anti-Chicken IgY (IgG) (H+L), highly cross-adsorbed (min X Bovine, Goat, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Sheep, and Syrian Hamster)
20275-1	50 uL	
20350	0.5 mL	Donkey Anti-Mouse IgG (H+L), highly cross-adsorbed (min X Bovine, Chicken, Goat, Guinea Pig, Horse, Human, Rabbit, Sheep, and Syrian Hamster)
20350-1	50 uL	
20351	0.5 mL	Donkey Anti-Rabbit IgG(H+L), highly cross-adsorbed (min X Bovine, Chicken, Goat, Guinea Pig, Horse, Human, Mouse, Rat, Sheep, and Syrian Hamster)
20351-1	50 uL	

Product Description

CF™350 antibodies are affinity-purified antibodies labeled with a blue fluorescent dye CF™350, one of an outstanding series of CF dyes developed by Biotium. CF™ dyes are an excellent alternative to Alexa Fluor® dyes, with superior brightness, photostability and specificity. CF™350 can be used as a direct replacement for Alexa Fluor® 350, AMCA, Cascade Blue®, Coumarin, DyLight™ 350.

Color and Form: Clear solution.

Concentration: 2 mg/mL in pH~7.4 PBS containing 50% glycerol, 2 mg/ml bovine serum albumin (IgG-free and protease-free) and 0.05% sodium azide.

Storage and Handling

Product is stable for at least 6 months at -20°C as an undiluted liquid. Storage of the antibody for more than a day at final working dilution is not recommended. Protect from light.


Recommended Dilution Range

1-10 µg/mL of the IgG conjugate for most applications (appropriate dilutions of the conjugate should be determined empirically). See other side for example staining protocols.

Spectral Properties

$\lambda_{abs}/\lambda_{em}$ = 347/448 nm in pH 7.4 PBS buffer (Figure 1).

Figure 1. Absorption/Emission Spectra of CF350 Conjugated Antibodies



Continued other side

General Protocols for Using CF™-labeled IgG Secondary Antibodies

Immunofluorescence Protocol for Microscopy

There are many methods for immunofluorescence staining. The protocol below is a general guideline for staining cells and should be optimized or modified to obtain the best results for each particular application.

1. Coverslip preparation for adherent cells

- 1.1 Culture cells on slide chambers or sterile glass coverslips (with poly-L-lysine coating if cells do not adhere well, see below). We recommend 18 x 18 mm square coverslips in 6-well plates or 4-well chamber slides.
- 1.2 Allow cells to adhere and treat as desired.
- 1.3 Rinse cells gently with PBS.

2. Coverslip preparation for non-adherent cells

- 2.1 Coat coverslips with 0.01% poly-L-lysine solution for 10 minutes at room temperature.
- 2.2 Aspirate the poly-L-lysine solution and allow coverslips to dry completely.
- 2.3 Centrifuge cells in medium and resuspend in PBS. Transfer cells to coverslips.
- 2.4 Incubate for 30-60 minutes. Check for adherence by microscope.

3. Fixation and Staining

- 3.1 Fix with 4% paraformaldehyde/PBS, 15 min.
- 3.2 Rinse twice with PBS to remove traces of fixative.
- 3.3 Permeabilize with 0.1 - 0.5% TritonX-100/PBS, 5-10 min.
- 3.4 Block with blocking agent such as with 5% BSA or normal goat serum in PBS, 30 min.
- 3.5 Dilute primary antibody in dilution buffer as recommended in the specific product's datasheet. Overlay enough diluted antibody to cover cells on coverslip (150-200 µL is usually sufficient to cover the surface area) or add to each chamber of the chamber slides. Keep slips covered or in a humidified chamber to avoid evaporation.
- 3.6 Rinse three times with PBS, 5 min each wash.
- 3.7 Dilute fluorescent secondary antibody in dilution buffer and incubate for 1 hour at room temperature. General range for IgG conjugates is between 1-10 µg/mL for most applications. Cell samples without primary antibody incubation is recommended for background control. Keep slips covered or in a humidified chamber to avoid evaporation.
- 3.8 Rinse three times with PBS, 5 min each wash.
- 3.9 Additional staining with fluorescent nuclear stains or phalloidins can be done at this step.
- 3.10 Invert each coverslip onto a pre-cleaned slide with fluorescence anti-fade mounting media. Seal edges with clear polish if desired.
- 3.11 Store slides in the dark at 4°C.

Staining Protocol for Flow Cytometry

There are many alternative procedures that can be used for specific staining experiments. The protocol below is a general guideline for flow cytometry and should be optimized or modified for each application.

- 1 Aliquot 1 X 10⁶ cells into 12 X 75 mm polypropylene tubes for flow cytometry.
- 2 For intracellular staining, cells can be fixed first to ensure stability of soluble antigens or antigens with short half-lives. We recommend a fix and perm kit from reliable manufacturers. Follow manufacturer's instructions.
- 3 Add the primary antibody or isotype control at the appropriate dilution to the assay tubes. Incubate according to manufacturer's instructions.
- 4 Rinse cells twice by centrifugation with 2-3 mL incubation buffer.
- 5 Decant supernatant and re-suspend the pellet in remaining volume of wash.
- 6 Add fluorescent secondary antibody and incubate for 20-30 minutes. General range for secondary antibodies is between 1-10 µg/mL for IgG conjugates for most applications.
- 7 Rinse cells twice by centrifugation with 2-3 mL incubation buffer. Centrifuge to collect cells after each wash. Decant supernatant.
- 8 Resuspend cells in 0.5 mL of diluent of choice to analyze on flow cytometer. Acquire data using the correct channel.

Tips and Hints:

- 1) No signal or weak fluorescence intensity may suggest the following: (a) insufficient antibody is present for detection, (b) intracellular target was not accessible, (c) excitation sources are not aligned, (d) target protein is not present or expressed at low levels, (e) fluorochrome has faded, and/or (f) primary and secondary antibodies are not compatible.
- 2) High fluorescence intensity may suggest the following: (a) antibody concentration is too high, (b) excess antibody was not washed away efficiently, and/or (c) blocking was inadequate. Increase antibody dilution and washes.

CF™-labeled antibodies can also be used for staining histological sections from paraffin-embedded or frozen tissues.

References

1. Donaldson, J.G. Immunofluorescence staining. (2001) Curr Protoc Cell Biol. Chapter 4: Unit 4.3.
2. Blose, S.H. and Feramisco, J.R. (1983) Fluorescent methods in the analysis of cell structure. Cold Spring Harbour Laboratory.

Useful websites:

www.chroma.com

Related Products

Please visit our website at www.biotium.com for more information about these related products

CF™ Dye Product	Application
Reactive dyes	Covalent conjugation of CF™ dyes to proteins, oligonucleotides and other biomolecules
Protein labeling kits	Covalent conjugation of CF™ dyes to proteins
Mix-n-Stain™ Antibody Labeling Kits	Single step covalent conjugation of 5-100 µg antibody to biotin or CF™ dyes in 30 minutes
Streptavidin Conjugates	Labeling biotin conjugates for microscopy, flow cytometry and Western blotting
Annexin V Conjugates	Labeling apoptotic cells for microscopy or flow cytometry
Phalloidin Conjugates	Labeling actin microfilaments for microscopy
α-Bungarotoxin	Labeling Ach-R for microscopy or flow cytometry
Lectin Conjugates	Labeling glycoproteins for microscopy or flow cytometry
EverBrite™ Mounting Medium	Antifade mounting medium available with or without DAPI
CoverGrip™ Coverslip Sealant	Nail polish alternative for sealing the edges of wet-mounted coverslips

Biotium products are for research use only. Not for use in diagnostic or therapeutic procedures. Our products are not available for resale or other commercial uses without a specific agreement from Biotium, Inc. We welcome inquiries about licensing the use of our dyes, trademarks or technologies. Please submit inquiries by e-mail to btinfo@biotium.com.

CF™ dye technology is covered by pending US and international patents. Alexa Fluor® and Cascade® Blue are registered trademarks of Invitrogen; DyLight® is a registered trademark of Thermo-Fisher Scientific.