



VWR® Bottle-Top Filtration Systems

Guaranteed Performance
with Maximum Throughput

Convenient, One-Time
Vacuum Connection

Complete System Allows
Filtration and Storage in
One Easy Step





VWR® Filtration Systems

VWR's new line of disposable, sterile laboratory filtration systems was designed for the preparation of buffers, tissue culture media, microbiological media, and other biological fluids. VWR's vacuum filtration systems feature large diameter Pall membranes for guaranteed performance and maximum throughput. The bottle-top filtration system design includes a receiver bottle with an ergonomic shape, molded finger grips, and a patented bottle cap. The centrifuge tube and funnel system allows the user to filter fluids directly into a 15mL or 50mL tube.

Each separate piece is compatible with the quick-connect vacuum pedestal, which allows for a convenient one-time vacuum hose connection, and one-handed operation. A cradle ring is also available for use with existing laboratory stands. Systems are available with 0.1µm (mycoplasma filtering grade), 0.2µm (sterilizing grade) or 0.45µm (clarification grade) asymmetric PES membranes. All items are made from FDA-grade, Class VI materials according to ISO 13485 standards and are gamma sterilized according to Vdmax 25 guidelines (SAL 10⁻⁶).

Filtrate Receiver/Storage Bottles

VWR's vacuum-rated, disposable round media bottles are available separately, for storing sterile solutions such as tissue culture media, serum, and buffers. Bottles feature a sloped design with a low center of gravity that improves stability on the work surface. The narrow bottle shoulder and molded finger grips enable the bottle to be gripped securely by gloved hands. Crystal clear bottles feature high-visibility metric graduations. Polypropylene caps feature a convenient gripping surface for easy handling and transporting. Caps allow for easy uncapping, tightening, and bottle stacking. The tops of caps also feature a flat surface for easy labeling.

Polycarbonate bottles can withstand up to three autoclave cycles and offer good alcohol resistance. They are stable up to 135°C (275°F). Polystyrene bottles are ideal for general purpose cell culture aqueous solutions, and are stable up to 80°C (176°F). PETG bottles are ideal for most general laboratory solutions, and offer gas impermeability as well as excellent impact resistance.






- 1 • Quick-Connect™ pedestal base or cradle ring allows for hands-free filtration
- 2 • Forward-facing, large, raised graduation marks allow for easy visualization of volumes
- 3 • Available in 0.1µm, 0.2µm, and 0.45µm PES asymmetric membrane filters
- 4 • Conveniently located vacuum on/off switch
- 5 • Available in a media storage bottle system or in a centrifuge funnel system
- 6 • Vacuum hose adapter allows for use without pedestal or cradle ring



Easy Initial Setup of VWR Filtration System


- Connect the vacuum hose just once to the back of the Quick-Connect pedestal
- Place filter on the Quick-Connect pedestal — the vacuum connects automatically
- Once fluid has been added to the funnel, turn the on/off switch to the ON position





VWR® Filtration Systems

- Feature large diameter Pall asymmetric PES membranes
- Quick-Connect™ pedestal allows for one-handed operation
- Complete system includes filtration funnel, media bottle or centrifuge tube. For Quick-Connect™ use pedestal, cradle ring, or hose adapter
- Sterile



Volume, mL	Cat. No.
0.1µm Centrifuge Tube Filtration Systems	
15	89220-714
50	89220-708
0.2µm Centrifuge Tube Filtration Systems	
15	89220-716
50	89220-710
0.45µm Centrifuge Tube Filtration Systems	
15	89220-718
50	89220-712
0.1µm Media Bottle Filtration Systems	
250	89220-694
500	89220-696
1000	89220-698
0.2µm Media Bottle Filtration Systems	
250	97066-200
500	97066-202
1000	97066-204
0.45µm Media Bottle Filtration Systems	
250	97066-206
500	97066-208
1000	97066-210



VWR® Filtration Funnels

- Compatible with the Quick-Connect™ pedestal
- Sterile



Volume, mL	Cat. No.
0.1µm Filtration Funnels	
15	89220-726
50	89220-720
250	89220-700
500	89220-702
1000	89220-704
0.2µm Filtration Funnels	
15	89220-728
50	89220-722
250	97066-212
500	97066-214
1000	97066-216
0.45µm Filtration Funnels	
15	89220-730
50	89220-724
250	97066-218
500	97066-220
1000	97066-222
Accessories	
Vacuum Filtration Cradle Ring	97066-234
Vacuum Filtration Pedestal	97066-232



VWR® Round Media Bottles

- Feature molded finger grips and a narrow shoulder for better gripping with gloved hands
- Unique polypropylene caps offer a flat writing surface and allow for bottles to be stacked
- Additional head space allows for 15% dilution without changing bottles

Capacity, mL (oz.)	Sterile	Cat. No.
PETG Bottles		
250 (8)	Yes	89166-098
250 (8)	No	89166-110
500 (16)	Yes	89166-100
500 (16)	No	89166-112
1000 (32)	Yes	89166-102
1000 (32)	No	89166-114
Polycarbonate Bottles		
250 (8)	Yes	89166-104
250 (8)	No	89166-122
500 (16)	Yes	89166-106

Capacity, mL (oz.)	Sterile	Cat. No.
Polycarbonate Bottles		
500 (16)	No	89166-124
1000 (32)	Yes	89166-108
1000 (32)	No	89166-126
Polystyrene Bottles		
250 (8)	Yes	89166-092
250 (8)	No	89166-116
500 (16)	Yes	89166-094
500 (16)	No	89166-118
1000 (32)	Yes	89166-096
1000 (32)	No	89166-120



1.800.932.5000 | vwr.com

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. © 2012 VWR International, LLC. All rights reserved.