

## PermaBlue/AP

(Alcohol & Xylene-substitute Compatible)

- Catalog No.:** K 051-110
- Intended Use:** Substrate/chromogen in conjunction with alkaline phosphatase based immunostaining or *in situ* hybridization systems.
- Introduction:** PermaBlue/AP is a substrate-chromogen system designed to be used for either IHC or ISH when utilizing alkaline phosphatase. PermaBlue/AP produces a vibrant blue color. It is insoluble in alcohol and xylene substitutes (both aliphatic hydrocarbon and citrus based); therefore sections can be dehydrated in alcohol, cleared in xylene substitute\*, and permanently mounted.
- Components:**
- (i) 110mL clear PermaBlue/AP Substrate Buffer.
  - (ii) 3mL concentrated PermaBlue/AP Chromogen.
  - (iii) One empty mixing dropping bottle.
- Storage:** Store at 2-8°C. Do not use beyond the expiration date stated on the label.
- Working Solution:** Aliquot 3mL of PermaBlue/AP Substrate Buffer in a mixing bottle. Add one drop (~20µl) of concentrated PermaBlue/AP Chromogen solution. Replace tip, mix, and allow solution to reach room temperature before using.  
*Note: The working PermaBlue/AP chromogen-substrate solution should be used within 25 minutes of preparation. Any solution not used during this period should be discarded.*
- Procedure:**
- i) After SA-alkaline phosphatase incubation, wash tissue sections with wash buffer.
  - ii) Wipe slides, removing excess buffer. Add enough drops of working PermaBlue/AP solution to cover tissue sections.
  - iii) Incubate for 5-25 minutes at room temperature. For optimal results, observe reaction under microscope for signal development. Once desired signal to noise ratio is achieved, stop reaction by rinsing slides with DI H<sub>2</sub>O.  
*NOTE: Increasing incubation temperature to 37°C will increase sensitivity and decrease needed incubation time.*
  - iv) Counter stain. Nuclear Fast Red provides good contrast. Wash with DI H<sub>2</sub>O.
  - v) Dehydrate sections in alcohol, clear in a xylene-substitute\* and mount in permanent mounting medium.

\*Notes: Use increasing concentrations of Ethanol up to 100% to dehydrate  
Use xylene-substitute instead of xylene

Alternatively, slides can be air dried (instead of dehydrated/cleared in alcohol and xylene-substitute). After rinsing off counter stain in DI H<sub>2</sub>O, leave slides on benchtop for at least 20 minutes to air dry. Then permanently mount.

### IVD: For In Vitro Diagnostic Use

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

**DBS**

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: [www.dbiosys.com](http://www.dbiosys.com) e-mail: [customersupport@dbiosys.com](mailto:customersupport@dbiosys.com)