

## Monoclonal Mouse Antibody to Human B Cell, CD20

<b>Catalog No.:</b>	PDM 004
<b>Intended Use:</b>	This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy. Clinical interpretation of staining results should be accompanied by histological studies with proper controls. Patients' clinical histories and other relevant diagnostic tests should be utilized by a qualified person(s) when evaluating and interpreting results.
<b>Immunogen:</b>	BALB/C mice were injected with human tonsil B cells.
<b>Clone:</b>	L26
<b>Isotype:</b>	IgG2a, kappa
<b>Format:</b>	This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded and acetone fixed cryostat tissue sections. No further titration is required.
<b>Staining Protocol:</b>	We suggest an incubation period of 30 minutes at room temperature. Optimal incubation conditions should be determined by the user based upon the fixation conditions and staining system employed. <u>Suitable for formalin fixed, paraffin embedded tissue sections.</u>
<b>Specificity:</b>	This antibody is specific to a 33 kD polypeptide present on the majority of B cells in peripheral blood and lymphoid tissue. It is also specific with a minor component of 30 kD. No reactivity with other hematopoietic cells has been observed.
<b>Positive Control:</b>	Tonsil
<b>Cellular Localization:</b>	Cell membrane
<b>Storage:</b>	Store at 2-8°C. Do not use beyond the expiration date stated on the label.
<b>References:</b>	i) Ishii et al. Clin Exp Immunol 58: 183, 1984. ii) Mason et al. Am J Pathol 136: 1215, 1990. iii) Cartun et al. Am J Pathol 129: 415, 1987.

### IVD: For In Vitro Diagnostic Use

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

**DBS**

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: [www.dbiosys.com](http://www.dbiosys.com) e-mail: [customersupport@dbiosys.com](mailto:customersupport@dbiosys.com)