

# Org Taqman QRT Kit


(Cat# M1187-100; One Step RNA Extraction from Cells and Taqman QRT PCR Kit; No Dye; Store at -20°C)

## I. Introduction:

Org Taqman QRT Kit offers a quick, simple and robust method to prepare template for performing Taqman Probe based quantitative real-time analysis directly from cultured cells, enabling reverse transcription of lysates from  $10^3$ - $10^5$  cultured cells without time-consuming and hazardous-chemicals-involved RNA extraction and purification steps. The kit includes reagents for cell lysis as well as gDNA removal. The lysis procedure simultaneously eliminates genomic DNA effectively in 12 min, without compromising RNA quality. The lysate can then be directly applied as template for One-Step QRT PCR, using the One Step Taqman QRT PCR Mastermix.

This One-Step Taqman QRT PCR system contains all necessary reagents for both reverse transcription and Taqman Probe based QPCR amplification to occur in a single qPCR reaction tube, including a QRT PCR Enzyme Mix in a proprietary buffer system to deliver precise and accurate sample analysis with high sensitivity and superb signal-to-noise ratio and an Taqman 2X QRT PCR Mastermix. The use of HotStart Taq DNA Polymerase in the enzyme mix significantly reduces non-specific PCR amplification. Our proprietary QRT PCR Enzyme Mix contains stabilizers and enhancers to optimize the two reactions in a real-time "single step". Coupled together, this complete system provides the ultimate convenience in generating consistent, reproducible, and accurate results from  $10^3$ - $10^5$  cells.

BioVision's Org Taqman QRT PCR Kit offers Extraction and preparing RNA templates directly from cultured cells to be applied as template for One-Step QRT PCR. Please refer to our QPCR Master Mix Selection Guide for selecting the appropriate QPCR formulation applicable to your particular instrument model.



## II. Application:

- Gene expression studies

## III. Package Contents (Org Taqman One Step QRT PCR Kit):

| Components                     | M1187-100 (25 lysis preps, 100 X 20 µl rxns) | Part Number |
|--------------------------------|--|-------------|
| Lysis Solution                 | 1.25 ml X 2                                  | M1187-XX-1  |
| Stop Solution                  | 300 µl                                       | M1187-XX-2  |
| Protease                       | 50 µl  | M1187-XX-3  |
| Protease Inhibitor             | 50 µl  | M1187-XX-4  |
| Taqman Master Mix-No Dye       | 1.25 ml                                      | M1187-XX-5  |
| QRT PCR Enzyme Mix (50X)       | 40 µl  | M1187-XX-6  |
| Nuclease-free H <sub>2</sub> O | 1 ml   | M1187-XX-7  |

## IV. User Supplied Reagents and Equipment:

- PCR Tubes
- Pipettes
- Water, Nuclease-free
- Primers (forward and reverse)
- Total RNA or poly(A) + mRNA

## V. Shipment and Storage:

Store all components at -20°C in a non-frost-free freezer. All components are stable for 1 year from the date of shipping when stored and handled properly. Avoid repeated freeze-thaw cycles to retain maximum performance. Briefly centrifuge small vials prior to opening.

## VI. Protocol:

- Thaw Lysis Solution and Stop Solution. Homogenize each solution gently but thoroughly.
- Prepare the following reactions for cell lysis:

| Components            | Volume |
|-----------------------|--------|
| $10^3$ - $10^5$ cells | 5 µl   |
| Protease              | 1 µl   |
| Lysis Solution        | 50 µl  |

Mix content by pipetting 35 µl of the mixture up and down 5 times and avoid creating bubbles. Incubate at 37°C for 10 min, then add the following to the tube:

|                    |      |
|--------------------|------|
| Protease Inhibitor | 1 µl |
| Stop Solution      | 5 µl |

Mix content by pipetting 35 µl of the mixture up and down 5 times. Incubate at room temperature for 2 min, then store the lysate on ice.

- The lysate is ready for QRT PCR setup. Prepare the following reaction mixture in a QPCR tube on ice:

| Components                 | Reaction Volume |       |       | Concentration |
|----------------------------|-----------------|-------|-------|---------------|
|                            | 10 µl           | 20 µl | 50 µl |               |
| Lysate from previous step* | 1 µl            | 2 µl  | 5 µl  | -             |

|                                     |             |             |             |            |
|-------------------------------------|-------------|-------------|-------------|------------|
| Taqman 2X QRT PCR Master Mix-No Dye | 5 µl        | 10 µl       | 25 µl       | 1X         |
| QRT PCR Enzyme Mix (50X)            | 0.2 µl      | 0.4 µl      | 1 µl        | 1X         |
| Forward Primer** (6 µM)             | 0.5 µl      | 1 µl        | 2.5 µl      | 300 nM     |
| Reverse Primer** (6 µM)             | 0.5 µl      | 1 µl        | 2.5 µl      | 300 nM     |
| Taqman Probe                        | Variable    | Variable    | Variable    | 100-300 nM |
| Nuclease-free H <sub>2</sub> O      | Up to 10 µl | Up to 20 µl | Up to 50 µl | -          |

Note:

1. As lysate is PCR-inhibitory in nature, the recommended volume of lysate used is 2%-10% of the reaction volume. If less than 10<sup>3</sup> cells are used to prepare the lysate, the QRT PCR system can tolerate more lysate, up to 20% of the reaction volume.
2. Gene specific primers must be used and amplicon should be
4. Gently mix and ensure all the components are at the bottom of the amplification tube. Centrifuge briefly if needed.
5. Program the thermal cycler so that cDNA synthesis is followed immediately by QPCR amplification.

| the thermal cycler so that cDNA synthesis is followed immediately by qPCR amplification. |  |          |           |
|--|--|----------|-----------|
| Steps  | Temperature                              | Duration | Cycle (s) |
| cDNA Synthesis   | 42°C                                     | 30 min   | 1         |
| Pre-Denaturation   | 95°C                                     | 10 min   | 1         |
| Denaturation   | 95°C                                     | 15 sec   | 40        |
| Annealing  | 60°C                                     | 60 sec   |           |
| Melt Curve   | According to the instrumental guidelines |          |           |

#### VII. General Notes:

- Minimize RNA degradation by keeping cells in PBS on ice before starting the cell lysis procedure.
- Do not vortex Stop Solution.
- Lysis Solution and Stop Solution must be at room temperature before proceeding to the lysis procedure. As cells settle quickly, thoroughly resuspend cells before withdrawing cell solution samples.
- (Optional) If setting up multiple reactions, prepare Protease/Lysis Solution premix for the number of reactions required, and then mix the premix solution with 5 µl of 10<sup>3</sup>-10<sup>5</sup> cells.
- (Optional) If setting up multiple reactions, prepare Protease Inhibitor / Stop Solution premix for the number of reactions required, and then mix the premix solution with the lysis reaction.
- As RNAs are poor templates for DNA polymerase, a Ct difference of 8-12 would be expected in QPCR between reactions containing RTase and those with no RTase.
- Lysates can be safely stored on ice for up to 1 hr after lysis. Alternatively, lysates can be stored at -80°C for a short period of time with a maximum of 1 freeze / thaw cycle. We highly recommend to use the lysates in downstream applications immediately after the 2 min termination.
- When handling 10<sup>3</sup> or more copy number of cells, ROX Referencing option could be turned off if applied by the QPCR instruments.

#### IX. Related Products:

| BV Product Name                        | BV Cat. No. |
|--|-------------|
| Two Step RT PCR Kits                   | M1160-M1161 |
| One Step RT PCR Kits                   | M1162-M1163 |
| First Strand cDNA Synthesis Kits       | M1164-M1167 |
| First Strand cDNA Synthesis Supermixes | M1167-M1169 |
| All-In-One RT Mastermixes              | M1170-M1172 |
| Reverse Transcriptases                 | M1173-M1174 |
| One Step Jade™ QRT PCR Kits            | M1175-M1182 |
| One Step Taqman QRT PCR Kits           | M1183-M1190 |

FOR RESEARCH USE ONLY! Not to be used on humans.