



## Product Specifications

# GENTAMYCIN SULFATE SOLUTION Tissue Culture Grade - 50mg/ml

<u>Catalog #</u>	<u>Size</u>
IB02030	20ml

### Physical Specifications

CAS#	N/A
Appearance	Clear colorless to light yellow solution
pH (@25°C)	4.5 ± 1.0
Specific Rotation	+107° - +121°

### Molecular Biology Specifications

Sterility	Negative
-----------	----------

### Recommended Use

Gentamycin Sulfate Solution is used as an inhibitor for the growth of a wide variety of gram-positive and gram negative microorganisms. For research use only.

### Storage

Store at +4°C. Warm to room temperature before opening. Keep tightly sealed.

### Warning

Not for drug use. Avoid contact with eyes, skin and clothing. Avoid breathing dust. Wash thoroughly with soap and water after handling. See Material Safety Data Sheet for additional information.