

**NBP1-80566****MEF2D Antibody****Unit Size**

0.1 ml

**Concentration**

0.2 mg/ml

**Storage**

Store at 4°C. Do not freeze.

**Applications**

WB, IP

**Species**

Human

**Host**

Rabbit

**Clonality**

Polyclonal

**Purity**

Immunogen affinity purified

**Preservative**

0.09% Sodium Azide

**Immunogen**

The epitope recognized by this antibody maps to a region between residue 425 and 475 of human Myocyte Enhancer Factor 2D using the numbering given in entry NP\_005911.1 (GeneID 4209).

**Recommended Dilutions**

Immunoprecipitation 2 - 10 ug/mg, Western Blot 1:2000 - 1:10000

**Buffer**


Tris-buffered Saline and 0.1% BSA

**Uses**

This antibody works in Western Blot, Immunoprecipitation.

**Images**

Western Blot: MEF2D Antibody [NBP1-80566] - Whole cell lysate from HeLa (15 and 50 mcg for WB; 1 mg for IP, 20% of IP loaded), 293T (T; 50 mcg) and Jurkat (J; 50 mcg) cells. Antibodies: Affinity purified rabbit anti-MEF2D antibody used for.... See more at [www.novusbio.com/NBP1-80566](http://www.novusbio.com/NBP1-80566)



Updated 3/9/2013 2.0

**Limitations:** This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee).**

**Novus USA**

8100 Southpark Way, A-8  
 Littleton, CO 80120  
 p) 888-506-6887  
 p) 303-730-1950  
 f) 303-730-1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

**Novus Canada**

461 North Service Road West,  
 Unit B37  
 Oakville, ON L6M 2V5  
 p) 855-668-8722  
 f) 905-827-6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

**Novus Europe**

12 Cambridge Science Park  
 Cambridge, UK CB4 0FQ  
 UK: p) +44-1223-426-001, f) +44-871-971-1635  
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79  
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340  
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31  
[europe@novusbio.com](mailto:europe@novusbio.com)

