

**NBP1-91170****TNF alpha ELISA Kit****Unit Size**

96 Tests

**Concentration**

Concentration is not relevant for this product.  
Please see the protocols for proper use of this product.

**Storage**

Store at 4C. Do not freeze.

**Applications**

ELISA

**Species**

Human

**Immunogen**

Suitable Sample Types: Cell culture supernatant, serum, plasma (EDTA, citrate, heparin)  
Sensitivity: 2.3 pg/mL Standard Curve Range: 7.8 - 500 pg/mL

**Recommended Dilutions**

ELISA

**Uses**

The human TNF-alpha ELISA is an enzyme-linked immunosorbent assay for the quantitative detection of human TNF-alpha

**Images**

This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: [www.novusbio.com/NBP1-91170](http://www.novusbio.com/NBP1-91170)

**Notes**

Components Included:

Aluminium pouch(es) with a Microwell Plate coated with monoclonal antibody to human TNF-alpha

Biotin-Conjugate anti-human TNF-alpha polyclonal antibody

Streptavidin-HRP

Human TNF-alpha Standard lyophilized, 1000 pg/mL upon reconstitution

Control high

, lyophilized

Control low

, lyophilized

Sample Diluent

Assay Buffer Concentrate ... See more at [www.novusbio.com/NBP1-91170](http://www.novusbio.com/NBP1-91170)

Updated 3/9/2013 2.0

**Limitations:** This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee).**

**Novus USA**

8100 Southpark Way, A-8  
Littleton, CO 80120  
p) 888-506-6887  
p) 303-730-1950  
f) 303-730-1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

**Novus Canada**

461 North Service Road West,  
Unit B37  
Oakville, ON L6M 2V5  
p) 855-668-8722  
f) 905-827-6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

**Novus Europe**

12 Cambridge Science Park  
Cambridge, UK CB4 0FQ  
UK: p) +44-1223-426-001, f) +44-871-971-1635  
DE: p) +49-6922-22340-60, f) +49-0800-58926-79  
IT: p) +39-02-4032-6786, f) +39-02-4032-6340  
FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31  
[europe@novusbio.com](mailto:europe@novusbio.com)

