

NBP1-27881**XPF Antibody (3F2/3)****Unit Size**

0.2 ml

Concentration

1.1 mg/ml

Storage

Store at -20C. Avoid freeze-thaw cycles.

Applications

WB, ELISA, ICC/IF

Species

Rat

Host

Mouse

Clonality

Monoclonal

Isotype

IgG1

Purity

Protein A purified

Preservative

No Preservative

Immunogen

A sequence to the C-terminus of XPF.

Recommended Dilutions

ELISA, Western Blot, Immunocytochemistry/Immunofluorescence

Buffer

PBS

Uses

ELISA, Immunofluorescence and Western Blot

References

Paul C, Nagano M, Robaire B. Aging Results in Differential Regulation of DNA Repair Pathways in Pachytene Spermatocytes in the Brown Norway Rat. Biol Reprod. 2011 Aug 24. [PMID: 21865553]

Images

This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: www.novusbio.com/NBP1-27881

Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

