

**NBP1-04354****SERPINA12 Antibody (1C4)****Unit Size**

0.1 ml

**Concentration**

1.0 mg/ml

**Storage**

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

**Applications**

WB, ELISA

**Species**

Human, Rat

**Host**

Mouse

**Clonality**

Monoclonal

**Isotype**

IgG1

**Purity**

Protein G purified

**Preservative**

0.1% Sodium Azide

**Specificity/Sensitivity**

Anti-human Vaspin mAb, clone 1C4, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human Vaspin protein.

**Immunogen**

Recombinant human Vaspin (21-414aa) purified from E. coli (NM\_173850).

**Recommended Dilutions**

ELISA, Western Blot 1:1000-1:3000

**Buffer**

PBS (pH 7.4)

**Images**

Western Blot: SERPINA12 Antibody (1C4) [NBP1-04354] Rat serum (15ug) was resolved by SDS-PAGE, transferred to NC membrane and probed with anti-human Vaspin (1:2000). Proteins were visualized using a goat anti-mouse secondary antibody conjugate... See more at [www.novusbio.com/NBP1-04354](http://www.novusbio.com/NBP1-04354)



Updated 3/9/2013 2.0

**Limitations:** This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee).**

**Novus USA**

8100 Southpark Way, A-8  
 Littleton, CO 80120  
 p) 888-506-6887  
 p) 303-730-1950  
 f) 303-730-1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

**Novus Canada**

461 North Service Road West,  
 Unit B37  
 Oakville, ON L6M 2V5  
 p) 855-668-8722  
 f) 905-827-6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

**Novus Europe**

12 Cambridge Science Park  
 Cambridge, UK CB4 0FQ  
 UK: p) +44-1223-426-001, f) +44-871-971-1635  
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79  
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340  
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31  
[europe@novusbio.com](mailto:europe@novusbio.com)

