

NBC1-18405**14-3-3 beta/alpha Protein****Unit Size**

0.1 mg

Concentration

1.0 mg/ml

Storage

Store at -80C. Avoid freeze-thaw cycles.

Applications

WB

Species

Human

Purity

Protein

Preservative

No Preservative

Specificity/Sensitivity

This protein is from amino acids 1-246

Immunogen

MTMDKSELVQKAKLAEQAERYDDMAAAMKAVTEQGHLSNEERNLLSVAYKNVVGARRS
 SWRVISSIEQKTERNEKKQQMGKEYREKIEAELQDICNDVLELLDKYLIPNATQPESKVFYK
 MKGDYFRYLSEVASGDNKQTTVSNSQAYQEAFAEISKEMQPTHPIRLGLALNFSVFYIEIL
 NSPEKACSLAKTAFDEAIAELDTLNEESYKDSTLIMQLLRDNLTLWTSENQGDGEGEGE
 N

Recommended Dilutions

Western Blot

Buffer


Liquid in 20mM Tris pH 8.0, 1mM EDTA, 50mM NaCl

Uses

Useful in Western Blot

Images

SDS-Page: 14-3-3 Beta Protein [NBC1-18405] - 14-3-3 beta, 28 kDa (246 aa), confirmed by MALDI-TOF with a purity of 90% by SDS - PAGE



Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

