

NB120-15307**metabotropic glutamate receptor 6/7 Antibody****Unit Size**

0.5 ml

Concentration

0.2 mg/ml

Storage

Store at 4C. Do not freeze.

Applications

IHC, IHC-P

Species

Human

Host

Rabbit

Clonality

Polyclonal

Purity

Immunogen affinity purified

Preservative

Sodium Azide

Immunogen

Synthetic peptide from the C-terminus of human mGluR6/7

Recommended Dilutions

Immunohistochemistry-Paraffin 1:100, Immunohistochemistry

Buffer


PBS containing 1% BSA pH 7.6

Uses

This antibody is useful in Immunohistochemical staining of formalin-fixed tissues.

Images

Immunohistochemistry-Paraffin: mGluR6/7 Antibody [NB120-15307] - Human brain stained with Anti-GluR6/7 antibody.



Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

