

NB110-8272**Staphylococcus Aureus Enterotoxin B Antibody (S222)****Unit Size**

0.2 mg

Concentration

Please see the vial label for concentration.

Storage

Store at 4C. Do not freeze.

Applications

IA, S-ELISA

Species

Bacteria

Host

Mouse

Clonality

Monoclonal

Isotype

IgG1

Purity

Protein G purified

Preservative

0.1% Sodium Azide

Specificity/Sensitivity

Clone S222 reacts with enterotoxin B of St. aureus but does not cross-react with St.aureus enterotoxins A, C1, C2, C3, D and E.

Immunogen

Recombinant St. aureus enterotoxin B expressed in E.coli

Recommended DilutionsImmunoassay
Sandwich ELISA**Buffer**

PBS with pH 7.4

UsesRecommended pair for immunoassay of St. aureus enterotoxin B (capture-detection):
NB110-8272 - [NB100-73025](#)**Images**This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: www.novusbio.com/NB110-8272

Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA
 8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com
Novus Canada
 461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com
Novus Europe
 12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com
