

NB100-74538

Somatostatin Receptor 3 Antibody

Unit Size

0.1 ml

Concentration

Concentrations vary lot to lot. See vial label for concentration.

Storage

Store at -20C. Avoid freeze-thaw cycles.

Applications

WB, ICC/IF, IHC, IHC-P

Species

Human

Host

Rabbit

Clonality

Polyclonal

Purity

Whole antisera

Preservative

0.05% Sodium Azide

Specificity/Sensitivity

Somatostatin Receptor 3

Immunogen

Synthetic peptide corresponding to residues C Q(a.a. 385) E R P P S R V A(a.a. 393) of human SSTR3.

Recommended Dilutions

Immunocytochemistry/Immunofluorescence 1:5000

Immunohistochemistry 1:2500

Immunohistochemistry-Paraffin 1:2500


Western Blot 1:5000

References

Arne G, Nilsson B, Dalmo J et al Gastrointestinal stromal tumors (GISTs) express somatostatin receptors and bind radiolabeled somatostatin analogs Acta Oncol 2012 Nov 1 [PMID: 23116418] (IHC, Human)

Images

Immunohistochemistry-Paraffin: Somatostatin Receptor 3 Antibody [NB100-74538] - Immunohistochemical staining of SSTR3 in human pancreatic islets using NB100-74538.



Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

