

NB100-74633**MARVELD2 Antibody****Unit Size**

0.1 mg

Concentration

1.0 mg/ml

Storage

Store at 4°C. Do not freeze.

Applications

WB, IP

Species

Human

Host

Rabbit

Clonality

Polyclonal

Purity

Immunogen affinity purified

Preservative

0.09% Sodium Azide

Immunogen

The immunogen recognized by this antibody maps to a region between residue 90 and 140 of human MARVEL domain containing 2 (Tricellulin) using the numbering given in entry NP_001033692.1 (GeneID 153562).

Recommended Dilutions

Immunoprecipitation 2-5 ug/mg lysate, Western Blot 1:2000-1:10000

Buffer

Tris-citrate/phosphate buffer, pH 7 to 8

Uses

This antibody is useful for western blot and immunoprecipitation.

References

Nguyen KH, Suzuki H, Wakasugi T et al Different expressions of erbB1/2 and tight junction proteins in hypertrophic inferior turbinates and nasal polyps Eur Arch Otorhinolaryngol 2012 Aug 28 [PMID: 22926991] (IHC-Fr, ICC/IF, Human)


Images

Western Blot: MARVELD2 Antibody [NB100-74633]

- Detection of Human MARVELD2/Tricellulin

HeLa whole cell lysate using NB100-74633.

MARVELD2/Tricellulin was also immunoprecipitated by rabbit anti-MARVELD2/Tricellulin antibody NB100-74632.



Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

