

**NB100-64510****Yellow Fever Virus Antibody (0G5)****Unit Size**

0.05 ml

**Concentration**

This product is unpurified. The exact concentration of antibody is not quantifiable.

**Storage**

Store at -20C. Avoid freeze-thaw cycles.

**Applications**

Func, ICC/IF

**Species**

Virus

**Host**

Mouse

**Clonality**

Monoclonal

**Isotype**

IgG2a

**Purity**

Ascites

**Preservative**

No Preservative

**Specificity/Sensitivity**

Specific for the envelope protein of the wild (Asibi) and vaccine strains of yellow fever virus. The clone has been reported to have neutralizing activity against the Asibi strain. No cross reactivity with other flaviviruses has been reported.

**Immunogen**

17D strain of yellow virus

**Recommended Dilutions**

Functional Assay, Immunocytochemistry/Immunofluorescence

**Uses**

Immunofluorescence and Functional Assays

**Images**

This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: [www.novusbio.com/NB100-64510](http://www.novusbio.com/NB100-64510)

Updated 3/9/2013 2.0

**Limitations:** This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee).**

**Novus USA**

8100 Southpark Way, A-8  
 Littleton, CO 80120  
 p) 888-506-6887  
 p) 303-730-1950  
 f) 303-730-1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

**Novus Canada**

461 North Service Road West,  
 Unit B37  
 Oakville, ON L6M 2V5  
 p) 855-668-8722  
 f) 905-827-6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

**Novus Europe**

12 Cambridge Science Park  
 Cambridge, UK CB4 0FQ  
 UK: p) +44-1223-426-001, f) +44-871-971-1635  
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79  
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340  
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31  
[europe@novusbio.com](mailto:europe@novusbio.com)

