

**NB100-62927**

## Staphylococcus Antibody

### Unit Size

0.5 ml

### Concentration

Please see the vial label for concentration.

### Storage

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

### Applications

ELISA

### Species

Bacteria

### Host

Rabbit

### Clonality

Polyclonal

### Purity

Affinity purified

### Preservative

0.09% Sodium Azide

### Specificity/Sensitivity

This product is specific for Staphylococcus enterotoxin B, an exotoxin excreted by Staphylococcus aureus that is the cause of classic food poisoning. It consists of 239 amino acids with a molecular weight of 28 kDa and is one of six antigenically distinct enterotoxin proteins that have been identified (A, B, C, D, E and G). This antibody does not react with other staphylococcal toxins.

### Immunogen

Purified toxin

### Recommended Dilutions

ELISA

### Buffer

Phosphate buffered saline

### Uses

This antibody is useful for ELISA

### Images

This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: [www.novusbio.com/NB100-62927](http://www.novusbio.com/NB100-62927)

Updated 3/9/2013 2.0

**Limitations:** This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee).**

#### Novus USA

8100 Southpark Way, A-8  
 Littleton, CO 80120  
 p) 888-506-6887  
 p) 303-730-1950  
 f) 303-730-1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

#### Novus Canada

461 North Service Road West,  
 Unit B37  
 Oakville, ON L6M 2V5  
 p) 855-668-8722  
 f) 905-827-6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

#### Novus Europe

12 Cambridge Science Park  
 Cambridge, UK CB4 0FQ  
 UK: p) +44-1223-426-001, f) +44-871-971-1635  
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79  
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340  
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31  
[europe@novusbio.com](mailto:europe@novusbio.com)

