

NB100-61104**Tryptophan Hydroxylase 1 Antibody****Unit Size**

0.1 ml

Concentration

Please see the vial label for concentration.

Storage

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

Applications

WB

Species

Rat

Host

Rabbit

Clonality

Polyclonal

Purity

Immunogen affinity purified

Preservative

No Preservative

Specificity/Sensitivity

Tryptophan Hydroxylase [Ser19]

Immunogen

Phosphopeptide corresponding to amino acid residues surrounding the phospho-Ser19 of rat tryptophan hydroxylase 2 (TPH2).

Recommended Dilutions

Western Blot 1:1000

Buffer

0.1 ml in 10mM HEPES [pH 7.5] containing 150 mM NaCl, 0.1 mg per ml BSA and 50% glycerol

Uses

This antibody is useful for Western Blot.

Images

Western Blot: Tryptophan Hydroxylase 1 [phospho Ser19] Antibody [NB100-61104] - Recombinant tryptophan hydroxylase incubated in the absence (Control) and presence of Ca²⁺/calmodulin dependent kinase II (CaMKII) showing specific immunolabelin... See more at www.novusbio.com/NB100-61104



Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

