

NB100-211**MKLP1 Antibody****Unit Size**

0.1 mg

Concentration

1.0 mg/ml

Storage

Store at 4C. Do not freeze.

Applications

WB, ICC/IF, IP

Species

Human

Host

Rabbit

Clonality

Polyclonal

Purity

Immunogen affinity purified

Preservative

0.09% Sodium Azide

Immunogen

The immunogen recognized by this antibody maps to a region between residues 1 and 50 of human Mitotic Kinesin-Like 1 using the numbering given in entry NP_004847.2 (GeneID 9493).

Recommended Dilutions

Immunocytochemistry/Immunofluorescence (PMID: 21652636), Immunoprecipitation 1-4mcg/mg of lysate, Western Blot 1:1000-1:10000

Buffer

Tris-citrate/phosphate buffer, pH 7 to 8

Uses


This antibody can be used for Immunofluorescence (PMID: 21652636).

References

Sinha S, Mondal G, Hwang EJ, et al. Von Hippel-Lindau gene product directs cytokinesis: a new tumor suppressor function. J Cell Sci;124(Pt 13):2132-2142. 2011 Jul 1. [PMID: 21470188]

Images

Western Blot: MKLP1 Antibody [NB100-211] - Whole cell lysate (10 ug/lane) from asynchronous (A) or mitotically-enriched (M) HeLa cells. Mitotically enriched HeLa cells were prepared by treatment with Nocodazole (500 ng/ml; 12h). Antibody us... See more at www.novusbio.com/NB100-211



Updated 3/9/2013 2.0

Limitations: This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit www.novusbio.com/guarantee.**

Novus USA

8100 Southpark Way, A-8
 Littleton, CO 80120
 p) 888-506-6887
 p) 303-730-1950
 f) 303-730-1966
novus@novusbio.com

Novus Canada

461 North Service Road West,
 Unit B37
 Oakville, ON L6M 2V5
 p) 855-668-8722
 f) 905-827-6402
canada@novusbio.com

Novus Europe

12 Cambridge Science Park
 Cambridge, UK CB4 0FQ
 UK: p) +44-1223-426-001, f) +44-871-971-1635
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com

