

**30920002****SM22 alpha Antibody****Unit Size**

0.05 mg

**Concentration**

Please see the vial label for concentration.

**Storage**

Store at -20C. Avoid freeze-thaw cycles.

**Applications**

WB, ELISA

**Species**

Human

**Host**

Rabbit

**Clonality**

Polyclonal

**Purity**

Affinity purified

**Preservative**

No Preservative

**Specificity/Sensitivity**

This product is specific for Human TAGLN.

**Immunogen**

This antibody was made against a protein fragment from the Middle Region

**Recommended Dilutions**

ELISA, Western Blot 1:5000-1:20000

**Buffer**

20 mM Potassium Phosphate, 150 mM Sodium Chloride, pH 7.0

**Uses**

This antibody is useful in ELISA and Western Blot.

**Images**

This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here: [www.novusbio.com/30920002](http://www.novusbio.com/30920002)

**Notes**

Manufactured by SDIX's proprietary Genomic Antibody Technology™. GAT [FAQs](#). These products are available in limited quantities. Please enquire about the availability of this product before purchasing.

Updated 3/9/2013 2.0

**Limitations:** This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. **For more information on our guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee).**

**Novus USA**

8100 Southpark Way, A-8  
 Littleton, CO 80120  
 p) 888-506-6887  
 p) 303-730-1950  
 f) 303-730-1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

**Novus Canada**

461 North Service Road West,  
 Unit B37  
 Oakville, ON L6M 2V5  
 p) 855-668-8722  
 f) 905-827-6402  
[canada@novusbio.com](mailto:canada@novusbio.com)

**Novus Europe**

12 Cambridge Science Park  
 Cambridge, UK CB4 0FQ  
 UK: p) +44-1223-426-001, f) +44-871-971-1635  
 DE: p) +49-6922-22340-60, f) +49-0800-58926-79  
 IT: p) +39-02-4032-6786, f) +39-02-4032-6340  
 FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31  
[europe@novusbio.com](mailto:europe@novusbio.com)

