



# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Revision Date 08/02/2012

Version 1.1

## SECTION 1. Identification

### Product identifier

Product number	109121
Product name	Potassium permanganate solution standardised against sodium thiosulfate $c(\text{KMnO}_4) = 0.02 \text{ mol/l}$ (0.1 N) TitriPUR®

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Reagent for analysis
-----------------	----------------------

### Details of the supplier of the safety data sheet

Company	EMD Millipore Corporation   290 Concord Road, Billerica, MA 01821, United States of America   SDS Phone Support: +1-978-715-1335   General Inquiries: +1-978-751-4321   Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)  e-mail: mm_sds@merckgroup.com
---------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Emergency telephone	800-424-9300 CHEMTREC (USA) +1-703-527-3887 CHEMTREC (International) 24 Hours/day; 7 Days/week
---------------------	------------------------------------------------------------------------------------------------------

## SECTION 2. Hazards identification

### GHS Classification

Chronic aquatic toxicity, Category 2, H411

For the full text of the H-Statements mentioned in this Section, see Section 16.

### GHS-Labeling

#### *Hazard pictograms*



#### *Hazard Statements*

H411 Toxic to aquatic life with long lasting effects.

#### *Precautionary Statements*

P273 Avoid release to the environment.

### OSHA Hazards

# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number	109121	Version 1.1
Product name	Potassium permanganate solution standardised against sodium thiosulfate $c(\text{KMnO}_4) = 0.02 \text{ mol/l}$ (0.1 N) TitriPUR®	

---

While this material is not considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200), this MSDS contains valuable information critical to the safe handling and proper use of the product. This MSDS should be retained and available for employees and other users of this product.

## Other hazards

None known.

---

## SECTION 3. Composition/information on ingredients

Chemical nature	Aqueous solution
-----------------	------------------

Remarks	No hazardous ingredients according to the OSHA Hazard Communication Standard 29CFR 19101200.
---------	----------------------------------------------------------------------------------------------

---

## SECTION 4. First aid measures

### Description of first-aid measures

#### *Inhalation*

After inhalation: fresh air.

#### *Skin contact*

After skin contact: wash off with plenty of water. Remove contaminated clothing.

#### *Eye contact*

After eye contact: rinse out with plenty of water.

#### *Ingestion*

After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell.

Never give anything by mouth to an unconscious person.

### Most important symptoms and effects, both acute and delayed

We have no description of any toxic symptoms.

### Indication of any immediate medical attention and special treatment needed

No information available.

---

## SECTION 5. Fire-fighting measures

### Extinguishing media

#### *Suitable extinguishing media*

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### Special hazards arising from the substance or mixture

Not combustible.

Ambient fire may liberate hazardous vapors.

### Advice for firefighters

# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number	109121	Version 1.1
Product name	Potassium permanganate solution standardised against sodium thiosulfate $c(\text{KMnO}_4) = 0.02 \text{ mol/l}$ (0.1 N) TitriPUR®	

---

## *Special protective equipment for fire-fighters*

In the event of fire, wear self-contained breathing apparatus.

## *Further information*

Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## SECTION 6. Accidental release measures

### Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel: Do not breathe vapors, aerosols. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

### Environmental precautions

Do not empty into drains.

### Methods and materials for containment and cleaning up

Cover drains. Collect, bind, and pump off spills.

Observe possible material restrictions (see sections 7 and 10).

Take up with liquid-absorbent material (e.g. Chemizorb®). Dispose of properly. Clean up affected area.

---

## SECTION 7. Handling and storage

### Precautions for safe handling

Observe label precautions.

### Conditions for safe storage, including any incompatibilities

Tightly closed.

Store at +15°C to +25°C (+59°F to +77°F).

---

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

Contains no substances with occupational exposure limit values.

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

### *Hygiene measures*

Change contaminated clothing. Wash hands after working with substance.

### *Eye/face protection*

Safety glasses

---

# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number	109121	Version 1.1
Product name	Potassium permanganate solution standardised against sodium thiosulfate $c(\text{KMnO}_4) = 0.02 \text{ mol/l}$ (0.1 N) TitriPUR®	

---

## *Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

Recommended:

full contact:

Glove material:	Nitrile rubber
Glove thickness:	0.11 mm
Break through time:	> 480 min

splash contact:

Glove material:	Nitrile rubber
Glove thickness:	0.11 mm
Break through time:	> 480 min

## *Respiratory protection*

required when vapors/aerosols are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

---

## SECTION 9. Physical and chemical properties

Physical state	liquid
Color	blue
Odor	odorless
Odor Threshold	No information available.
pH	ca. 8 at 68 °F ( 20 °C)
Melting point	No information available.
Boiling point	No information available.
Flash point	No information available.
Evaporation rate	No information available.
Flammability (solid, gas)	not applicable
Lower explosion limit	No information available.
Upper explosion limit	No information available.
Vapor pressure	No information available.

## MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number	109121	Version 1.1
Product name	Potassium permanganate solution standardised against sodium thiosulfate $c(\text{KMnO}_4) = 0.02 \text{ mol/l}$ (0.1 N) TitriPUR®	

Relative vapor density	No information available.
Relative density	$1.01 \text{ g/cm}^3$ at 68 °F ( 20 °C)
Water solubility	at 68 °F ( 20 °C) soluble
Partition coefficient: n-octanol/water	No information available.
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	No information available.

### SECTION 10. Stability and reactivity

#### Reactivity

See below

#### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

#### Possibility of hazardous reactions

The generally known reaction partners of water.

#### Conditions to avoid

no information available

#### Incompatible materials

no information available

#### Hazardous decomposition products

no information available

### SECTION 11. Toxicological information

#### Information on toxicological effects

##### *Likely route of exposure*

Eye contact, Skin contact

##### *Specific target organ systemic toxicity - single exposure*

The substance or mixture is not classified as specific target organ toxicant, single exposure.

##### *Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

##### *Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

#### Carcinogenicity

## MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number	109121	Version 1.1
Product name	Potassium permanganate solution standardised against sodium thiosulfate $c(\text{KMnO}_4) = 0.02 \text{ mol/l}$ (0.1 N) TitriPUR®	

---

IARC	No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
OSHA	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
ACGIH	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

### Further information

Quantitative data on the toxicity of this product are not available.

Further toxicological data:

Hazardous properties cannot be excluded, but are relatively improbable due to the low concentration of the dissolved substance(s).

Further data:

Handle in accordance with good industrial hygiene and safety practice.

---

## SECTION 12. Ecological information

### Ecotoxicity

No information available.

### Persistence and degradability

No information available.

### Bioaccumulative potential

No information available.

### Mobility in soil

No information available.

### Other adverse effects

*Additional ecological information*

Biological effects:

Bactericidal effect.

Further information on ecology

Discharge into the environment must be avoided.

# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number	109121	Version 1.1
Product name	Potassium permanganate solution standardised against sodium thiosulfate $c(\text{KMnO}_4) = 0.02 \text{ mol/l}$ (0.1 N) TitriPUR®	

---

---

## SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

---

## SECTION 14. Transport information

### Land transport (DOT)

UN number	UN 3082
Proper shipping name	ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S. ( POTASSIUM PERMANGANATE)
Class	9
Packing group	III
Environmentally hazardous	--

### Air transport (IATA)

UN number	UN 3082
Proper shipping name	ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S. ( POTASSIUM PERMANGANATE SOLUTION)
Class	9
Packing group	III
Environmentally hazardous	--
Special precautions for user	no

### Sea transport (IMDG)

# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)

Product number	109121	Version 1.1
Product name	Potassium permanganate solution standardised against sodium thiosulfate c(KMnO <sub>4</sub> ) = 0.02 mol/l (0.1 N) TitriPUR®	

UN number	UN 3082
Proper shipping name	ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S. ( POTASSIUM PERMANGANATE SOLUTION)
Class	9
Packing group	III
Environmentally hazardous	--
Special precautions for user	yes
EmS	F-A S-F

## SECTION 15. Regulatory information

### United States of America

#### OSHA Hazards

No OSHA Hazards

This information is based on 29 CFR 1910.1200 criteria prior to adoption of the GHS, and may deviate from the GHS information on the label and in section 2.

#### SARA 311/312 Hazards

No SARA Hazards

### US State Regulations

#### Massachusetts Right To Know

Remarks

No components are subject to the Massachusetts Right to Know Act.

#### Pennsylvania Right To Know

*Ingredients*

water

#### New Jersey Right To Know

*Ingredients*

water

#### California Prop 65 Components

This product does not contain any chemicals known to the State of California to cause cancer, birth, or any other reproductive defects.

#### Notification status



# MATERIAL SAFETY DATA SHEET

according to the Hazard Communication Standard (29 CFR 1910.1200)


Product number	109121	Version 1.1
Product name	Potassium permanganate solution standardised against sodium thiosulfate $c(\text{KMnO}_4) = 0.02 \text{ mol/l}$ (0.1 N) TitriPUR®	

TSCA: On TSCA Inventory

DSL: All components of this product are on the Canadian DSL list.

## SECTION 16. Other information

### National Fire Protection Association (U.S.A)



### Training advice

Provide adequate information, instruction and training for operators.

### Full text of H-Statements referred to under sections 2 and 3.

H411 Toxic to aquatic life with long lasting effects.

### Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*