

Datasheet

ALOX12 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00000239-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ALOX12 protein.

Immunogen: ALOX12 (NP_000688, 1 a.a. ~ 663 a.a) full-length human protein.

Sequence:

MGRYRIRVATGAWLFSGSYNRVQLWLVGTRGEAELE
LQLRPARGEEEEFDHDAEDLGLLQFVRLRKHHWLVD
DAWFCDRITVQGPACAEVAFPCYRWVQGEDILSLPE
GTARLPGDNALDMFQKHREKELKDRQQIYCWTWKE
GLPLTIAADRKDDLPPNMRFHEEKRLDFEWTLKAGAL
EMALKRVYTLSSWNCLEDFDQIFWGQKSALAEKVRQ
CWQDDELFSYQFLNGANPMLLRSTSLPSRLVLPSPG
MEELRAQLEKELQNGSLFEADFILLDGIPANVIRGEKQ
YLAAPLVMLKMEPNGKLQPMVIQPPNPSSPTPTLFL
PSDPPLAWLLAKSWVRNSDFQLHEIQYHLLNTHLVAE
VIAVATMRCLPGLHPIFKFLPHIRYTMETRARTQLIS
DGGIFDKAVSTGGGGHVQLLRRAAQLTYCSLCPDD
LADRGLLGLPGALYAHDAALRLWEIARYVEGIVHLFYQR
DDIVKGDPELQAWCREITEVGLCQAQDRGFPVSFQSQ
SQLCHFLTMCVFTCTAQHAAINQGQLDWYAWVPNAP
CTMRMPPTTKEDVTMATVMGSLPDVRQACLQMAIS
WHLRRQPDMPVPLGHHKEKYFSGPKPKAVLNQFRD
LEKLEKEITARNEQLDWPYEYLPSCIENSVTI

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 239

Gene Symbol: ALOX12

Gene Alias: 12-LOX, 12S-LOX, LOG12