

Datasheet

ALPPL2 polyclonal antibody (A01)

Catalog Number: H00000251-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant ALPPL2.

Immunogen: ALPPL2 (NP_112603, 365 a.a. ~ 454 a.a) partial recombinant protein with GST tag.

Sequence:

SEEDTLSLVTADHSHVFSFGGYPLRGSSIFGLAPGKAR
DRKAYTVLLYGNGPGYVLKDGARPDVTESESGSPEYR
QQSAVPLDGETHAGE

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 251

Gene Symbol: ALPPL2

Gene Alias: ALPG, ALPPL, GCAP

Gene Summary: There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The product of this gene is a membrane bound glycosylated enzyme, localized to testis, thymus and certain germ cell tumors, that is closely related to both the placental and intestinal forms of alkaline phosphatase. [provided by RefSeq]