

Datasheet

UNC5D polyclonal antibody (A01)

Catalog Number: H00137970-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant UNC5D.

Immunogen: UNC5D (NP_543148, 31 a.a. ~ 130 a.a) partial recombinant protein with GST tag.

Sequence:

TDNGEALPESIPSAPGTLPHFIEEPDDAYIIKSNPIALRC
KARPAMQIFFKCNGEVWHQNEHVSEETLDESSGLKV
REVFINVTRQQVEDFHGPEDYWCQ

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 137970

Gene Symbol: UNC5D

Gene Alias: FLJ16019, KIAA1777, PRO34692, Unc5h4