

Datasheet

PASD1 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00139135-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human PASD1 protein.

Immunogen: PASD1 (AAH40301.1, 1 a.a. ~ 773 a.a) full-length human protein.

Sequence:

MKMRGEKRRDKVNP KSSQRKLNWIP SFPTYDYFNQV
TLQLLDGFMITLSTDGVIICVAENISSLLGHLPAEIVGKK
LLSLLPDEEKDEVYQKIILKFLLNSETHIEFCCHLKRG
VEHGDSSAYENVKFIVNVRDICNEFPVVSGLFSSHL
ADFAACVPQEDRLYL VGNVCILRTQLLQQLYTSKAVSD
EAVLTQDSDEEPFVGELSSSQGQRGHTSMKAVYVEP
AAAAAAAAISDDQIDIAEVEQYGPQENVHMFVSDSTY
CSSTVFLDTMPESPALS LQDFRGEPEVNPLYRADPVD
LEFSVDQVDSVDQEGPMDQQDPENPVAPLDQAGLMD
PVPDPSVDLGAAGASAQPLQPSSPVAYDIISQELELM
KKLKEQLEERTWLLHDAIQNQNALELMDHLQKQPN
TLRHVVIPDLQSSEAVPKKQKQHAGQVKRPLPHPKD
VKGFCGLSLSNSLKN TGELQEPCVAFNQQLVQQEQ
HLKEQQRQLREQLQQLREQRKVQKQKKMQEKKLQ
QKMQEKKLQEQRQK KKKLQERKKWQGQMLQKEP
EEEQKQQLQEQLKHN VIVGNERNVQICLQNPRDVS
PLCNHPVRFLQAQPIVPVQRAAEQQPSGFYQDENG
QQEDESQSFYPEAYQGPPVNQLPLIDTSNSEAISSSI
PQFPITSDSTISTLETPQDYIRLWQELSDSLGPVVQVNT
WSCDEQGT LHGQPTYHQVQVSEVGVEGPPDPQAFQ
GPAAYQPDQMRSAEQTRLMPAEQRDSNKPC

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 139135

Gene Symbol: PASD1

Gene Alias: -