

Datasheet

RP11-311P8.3 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00139596-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human RP11-311P8.3 protein.

Immunogen: RP11-311P8.3 (NP_659489, 1 a.a. ~ 309 a.a) full-length human protein.

Sequence:

MATELQCPDSMPCHNQVNSASTPSPEQLRPGDLILD
HAGGNRASRAKVILLTGYAHSSSLPAELDSGACGGSSL
NSEGNSGSGDSSSYDAPAGNSFLEDCELSRQIGAQLK
LLPMNDQIRELQTIIRDKTASRGDFMFSADRLIRLVVEE
GLNQLPYKECMVTTPTGYKYEGVKFEKGNGVVSIMRS
GEAMEQGLRDCCRSIRIGKILIQSDEETQRAKVYYAKF
PPDIYRRKVLLMYPILSTGNTVIEAVKVLIEHGVQPSVIL
LSLFSTPHGAKSIIQEFPEITILTTEVHPVAPTHFGQKYF
GTD

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 139596

Gene Symbol: UPRT

Gene Alias: DKFZp781E1243, FUR1, MGC23937, UPP