

## Datasheet

### ASB8 polyclonal antibody (A01)

**Catalog Number:** H00140461-A01

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse polyclonal antibody raised against a partial recombinant ASB8.

**Immunogen:** ASB8 (NP\_077000, 179 a.a. ~ 288 a.a)  
partial recombinant protein with GST tag.

**Sequence:**

VINLIGQTPISRLVALLVRGLGTEKEDSCFELLHRAVGH  
FELRKNGTMPREVARDPQLCEKLTVLC SAPGTLKTLA  
RYAVRRSLGLQYLPDAVKGLPLPASLKEYLLLLLE

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** 50 % glycerol

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 140461

**Gene Symbol:** ASB8

**Gene Alias:** FLJ21255, MGC5540, PP14212