

Datasheet

ASB8 monoclonal antibody (M05), clone 1H1

Catalog Number: H00140461-M05

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant ASB8.

Clone Name: 1H1

Immunogen: ASB8 (NP_077000, 179 a.a. ~ 288 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

VINLIGQTPISRLVALLVRGLGTEKEDSCFELLHRAVGH
FELRKNGTMPREVARDPQLCEKLTVLCSAPGTLKTLA
RYAVRRSLGLQYLPDAVKGLPLPASLKEYLLLLLE

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 140461

Gene Symbol: ASB8

Gene Alias: FLJ21255, MGC5540, PP14212