

Datasheet

ZFP28 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00140612-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZFP28 protein.

Immunogen: ZFP28 (AAH99903.1, 1 a.a. ~ 349 a.a)
full-length human protein.

Sequence:

MRGAASASVHEPTPLPGRGAPRTKPRAGRGPTVGTP
ATLALPARGRPRSRNGLASKGQRGAAPTGPGRALP
SRDTALPQERNKKLEAVGTGIEPKAMSQGLVTFGDVA
VDFSQEEWEWLNPIQRNLYRKVMLENYRNLASLGCV
SKPDVISSLEQGKEPWTVKRKMTRAWCPDLKAVWKIK
ELPLKKDFCEGKLSQAVITERLTSYNLEYSLLGEHWDY
DALFETQPGLVTIKNLAVDFRQQLHPAQKNFCKNGIW
ENNSDLGSAGHCVAKPDVLSLLEQEKEPWMVKRELT
GSLFSGECRRRAWCGKSTARESGAFTKKAASLIYSRK
LYPGSLGLERKLNSGSSK

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 140612

Gene Symbol: ZFP28

Gene Alias: DKFZp686K03205, KIAA1431,
MGC104540, mkr5