

Datasheet

ZBTB46 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00140685-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZBTB46 protein.

Immunogen: ZBTB46 (AAH52269.1, 1 a.a. ~ 589 a.a) full-length human protein.

Sequence:

MNNRKEDMEIASHYRHLLRELNEQRQHGVLCDVCVV
VEGKVFKAHKNVLLGSSRYFKTLYCQVQKTSEQATVT
HLDIVTAQGFKAIIDFMYSAPHLALTSRNVIEVMSAASFL
QMTDIVQACHDFIKAALDISIKSDASDELAEFEIGASSS
SSTEALISAVMAGRSISPWLARRTSPANSSGDSAIASC
HDGGSSYGKEDQEPKADGPDDVSSQPLWPGDVGYG
PLRIKEEQVSPSQYGGSELPKADGAVQNSFSEQSAG
DAWQPTGRRKNRKNKETVRHITQQVEDDSRASSPVP
SFLPTSGWPFSSRDSNADLSVTEASSSDSRGERAELY
AQVEEGLLGGEASYLGPPLTPEKDDALHQATAVANLR
AALMSKNSLLSLKADVLGDDGSLLEFYLPRGAHSLSLN
EFTVIRKKFKCPYCSFSAMHQCILKRHMRSHTGERPY
PCEICGKKFTRREHMKRHTLVHSDKKYVCKVCSRVF
MSAASVGIRHGSRRHGVCTDCAGRGMAGPLDHGGG
GGEGSPEALFPGDGPYLEDPEDPRGEAEELGEDDEG
LAPEDALLADDKDEEDSPRPRSPPGPDKDFAWLS

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 140685

Gene Symbol: ZBTB46

Gene Alias: BTBD4, FLJ13502, RINZF, ZNF340, dJ583P15.7, dJ583P15.8