

## Datasheet

### RAB40A purified MaxPab mouse polyclonal antibody (B01P)

**Catalog Number:** H00142684-B01P

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse polyclonal antibody raised against a full-length human RAB40A protein.

**Immunogen:** RAB40A (NP\_543155.2, 1 a.a. ~ 277 a.a) full-length human protein.

**Sequence:**

MSAPGSPDQAYDFLLKFLLVGDRDVGKSEILESQDG  
AAESPYSHLGGIDYKTTTILLDGQRVKLKLWDTSGQGR  
FCTIFRSYRGAQGVILVYDIANRWSFEGMDRWIKKIE  
EHAPGVPKILVGNRLHLAFKRQVPREQAQAYAERLGV  
TFFEVSPLCNFNIESFTELARIVLLRHRMNWLGRPSKV  
LSLQDLCCRTIVSCTPVHLVDKLPLPSTLRSHLKSFSM  
AKGLNARMMRGLSYSLTTSSTHKSSLCKVEIVCPPQS  
PPKNCTRNSCKIS

**Host:** Mouse

**Reactivity:** Human

**Applications:** WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 142684

**Gene Symbol:** RAB40A

**Gene Alias:** MGC142061, RAR2A, Rar-2