

Datasheet

FLJ38725 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00144811-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human FLJ38725 protein.

Immunogen: FLJ38725 (AAH35749.1, 1 a.a. ~ 430 a.a) full-length human protein.

Sequence:

MAEAVLIDLFGLKLNSQKNCHQTLLKTLNAVQYHHA
AKFLCIMCCSNISYERDGEQDNCEIETSNGLSALLEEF
EIVSCPSMAATLYTIKQKIDKLNLSIKVIVPRHRKTLMK
AFIDQLFTDVYNFEFEDLQVTFRGGLFKQSIEINVITAQ
ELRGIQNEIETFLRSLPALRGKLTITSSLIPDIFIHGFTTR
TGGISYIPTLSSFNLFSSSKRRDPKVVVQENLRRLANA
AGFNVEKFYRIKTHHSNDIWIMGRKEPDSYDGITTNR
GVTIAALGADCIPVFADPVKKACGVAHAGWKGTLLGV
AMATVNMAMIAEYGCSELDIVVVLGPSVGPCCFITLPRES
AEAFHNLHPACVQLFDSPNPCIDIRKATRILLEQGGILP
QNIQDQNQDLNLCTSCHPDKFFSHVRDGLNFGTQIGFI
SIKE

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 144811

Gene Symbol: C13orf31

Gene Alias: DKFZp686D11119, FLJ38725