

Datasheet

GSC polyclonal antibody (A01)

Catalog Number: H00145258-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant GSC.

Immunogen: GSC (NP_776248, 151 a.a. ~ 257 a.a)
partial recombinant protein with GST tag.

Sequence:

LLNQLHCRRKRRHRTIFTDEQLEALENLFQETKYPDVG
TREQLARKVHLREEKVEVWFKNRRRAKWRRQKRSSSE
ESENAEKWNKTSSSKASPEKREEEGKSDLSDS

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 145258

Gene Symbol: GSC

Gene Alias: -

Gene Summary: This gene encodes a member of the bicoid subfamily of the paired (PRD) homeobox family of proteins. The encoded protein acts as a transcription factor and may be autoregulatory. A similar protein in mice plays a role in craniofacial and rib cage development during embryogenesis. [provided by RefSeq]