

Datasheet

MGC5987 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00145553-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MGC5987 protein.

Immunogen: MGC5987 (NP_612485.2, 1 a.a. ~ 176 a.a) full-length human protein.

Sequence:

MARLPKLAVFDLDYTLWPFWVDTHVDPPFHKSSDGTV
RDRRGQDVRLYPEVPEVLKRLQSLGVPGAAASRTSEI
EGANQLLELFDLFRYFVHREIYPGSKITHFERLQQKTGI
PFSQMIFFDDERRNIVDVSKLGVTCIHQNGMNLQTLS
QGLETFAKAQTGPLRSSLEESPFEA

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 145553

Gene Symbol: MDP-1

Gene Alias: MGC5987