

Datasheet

HAPLN3 polyclonal antibody (A01)

acid binding and cell adhesion. [provided by RefSeq]

Catalog Number: H00145864-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length recombinant HAPLN3.

Immunogen: HAPLN3 (AAH62320.1, 18 a.a. ~ 360 a.a) full-length recombinant protein with GST tag.

Sequence:

LPFYNGFYYSNSANDQNLGNHGHGKDLLNGVKLVVETP
EETLFTYQGASVILPCRYRYEPALVSPRRVRVKWWKL
SENGAPEKDVLAIGLRHRSFGDYQGRVHLRQDKEH
DVSLEIQDLRLLEDYGRYRCEVIDGLEDESGLVELELRG
VVFPYQSPNGRYQNFHEGQQVCAEQAAVASFEQL
FRAWEEGLDWCNAGWLQDATVQYPIMLPRQPCGGP
GLAPGVRSYGPRHRLHRYDVFCFATALKGRVYYLEH
PEKLTLEAREACQEDDATIAKVGQLFAAWKFHGLDR
CDAGWLADGSVRYPVVHPHPNCGPPEPGVRSFGFPD
PQSRLYGVCYRQH

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 145864

Gene Symbol: HAPLN3

Gene Alias: EXLD1, HsT19883

Gene Summary: This gene belongs to the hyaluronan and proteoglycan binding link protein gene family. The protein encoded by this gene may function in hyaluronic