

Datasheet

HAPLN3 purified MaxPab mouse polyclonal antibody (B01P)

protein encoded by this gene may function in hyaluronic acid binding and cell adhesion. [provided by RefSeq]

Catalog Number: H00145864-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human HAPLN3 protein.

Immunogen: HAPLN3 (NP_839946.1, 1 a.a. ~ 360 a.a) full-length human protein.

Sequence:

MGLLLLVPLLLLPGSYGLPFYNGFYYSNSANDQNLGN
GHGKDLLNGVKLVVETPEETLFTYQGASVILPCRYRYE
PALVSPRRVRVKWWKLENGAPEKDVLAIGLRHRSF
GDYQGRVHLRQDKEHDVSLEIQDLRLDYGRYRCEVI
DGLEDESGLVELELRGVVFPYQSPNGRYQFNFHEGQ
QVCAEQAAVVASFEQLFRAWEEGLDWCNAGWLQDA
TVQYPIMLPRQPCGGPGLAPGVRSYGPRHRRRLHRYD
VFCFATALKGRVYYLEHPEKLTLEAREACQEDDATIA
KVGQLFAAWKFFHGLDRCDAGWLADGSVRYPVHHPHP
NCGPPEPGVRSFGFDPQSRLYGVCYRQH

Host: Mouse

Reactivity: Human

Applications: WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 145864

Gene Symbol: HAPLN3

Gene Alias: EXLD1, HsT19883

Gene Summary: This gene belongs to the hyaluronan and proteoglycan binding link protein gene family. The