

Datasheet

MGC23280 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00147015-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MGC23280 protein.

Immunogen: MGC23280 (NP_653284, 1 a.a. ~ 327 a.a) full-length human protein.

Sequence:

MTALELARRGARVVLACRSQERGEAAAFDLRQESGN
NEVIFMALDLASLASVRAFATAFLSSEPRLDILIHNAGIS
SCGRTREAFNLLLRVNHIGPFLLTHLLLPCLKACAPSR
VVVVASAAHCRGRDLDFKRLDRPVVGWRQELRAYADT
KLANVLFARELANQLEATGVTCYAAHPGPVNSELFRL
HVPGWLRPLLRLPLAWLVLRAPRGGAQTPLYCALQEGI
EPLSGRYFANCHVEEVPPAARDDDRAAHLWEASKRL
AGLGPGEDAEPDEDPQSEDSEAPSSLSTPHPEEPTVS
QPYPSPQSSPDLSKMTHRIQAKVEPEIQLS

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 147015

Gene Symbol: DHRS13

Gene Alias: MGC23280, SDR7C5