

Datasheet

APCDD1 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00147495-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human APCDD1 protein.

Immunogen: APCDD1 (NP_694545.1, 1 a.a. ~ 514 a.a) full-length human protein.

Sequence:

MSWPRRLRLRYLFPALLHGLGEGSALLHPDSRSHPR
SLEKSAWRAFKEQCHHMLKHLHNGARITVQMPPTIE
GHWVSTGCEVRSGPEFITSYRFYHNNTFKAYQFYYG
SNRCTNPTYTLIRGKIRLRQASWIIRGGTEADYQLHNV
QVICHTEAVAELGQQVNRTPGFLADGGPWVQDVA
YDLWREENGCECTKAVNFAMHELQLIRVEKQYLHNL
DHLVEELFLGDIHTDATQRMFYRPSSYQPPLQNAKNH
DHACIACRIIYRSDEHPPILPPKADLTIGLHGEWVSQR
CEVRPEVFLTRHFIFHDNNNTWEGHYHYSDPVCKH
PTFSIYARGRYSRGVLSSRVMMGGTEFVFKVNHMKVTP
MDAATASLLNVFNGNECGAEGSWQVGIQQDVHTNG
CVALGIKLPHTYEYIFKMEQDARGRYLLFNGQRPSDG
SSPDRPEKRATSYQMPLVQCASSPRAEDLAEDSGS
SLYGRAPGRHTWSLLLAALACLVPLLHWNIRR

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 147495

Gene Symbol: APCDD1

Gene Alias: B7323, DRAPC1, FP7019

References:

1. APCDD1 is a novel Wnt inhibitor mutated in hereditary hypotrichosis simplex. Shimomura Y, Agalliu D, Vonica A, Luria V, Wajid M, Baumer A, Belli S, Petukhova L, Schinzel A, Brivanlou AH, Barres BA, Christiano AM. Nature. 2010 Apr 15;464(7291):1043-7.