

Datasheet

HIPK4 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00147746-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human HIPK4 protein.

Immunogen: HIPK4 (NP_653286.2, 1 a.a. ~ 616 a.a) full-length human protein.

Sequence:

MSTIQSETDCYDIEVLGKGTGGEVAKGWRRSTGEMV
AIKILKNDAYRNRIIKNELKLLHCMRGLDPEEAHVIRFLE
FFHDALKFYLVFELLEQNLFEFQKENNFAPLPARHIRT
TLQVLTALARLKLAIHADLKPENIMLVDQTRCPFRVK
VIDFGSASIFSEVRYVKEPYIQSRFYRAPEILLGLPFCEK
VDVWSLGCVMAELHLGWPLYPGNNEYDQVRYICETQ
GLPKPHLLHAACKAHFFKRNPHPDAAANPWQLKSSA
DYLAETKVRPLERRKYMLKSLDQIETVNGGSVASRLTF
PDREALAEHADLKSMVELIKRMLTWESHERISPSAALR
HPFVSMQQLRSAHETTHYYQLSLRSYRLSLQVEGKPP
TPVVAEDGTPYYCLAEKEAAGMGSVAGSSPFFREE
KAPGMQRAIDQLDDLSLQEAGHGLWGETCTNAVSDM
MVPLKAAITGHHVPDSGPEPILAFYSSRLAGRHKARKP
PAGSKSDSNFSNLIRLSQVSPEDDRPCRGSSWEEGE
HLGASAEPLAILQRDEDGPNIDNMTMEAERPDPFLFD
PSSCPGEWLSEPDCTLESVRGPRAQGLPPRRSHQHGH
PPRGATSFLQHVTGHH

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 147746

Gene Symbol: HIPK4

Gene Alias: FLJ32818