

## Datasheet

### HIPK4 purified MaxPab rabbit polyclonal antibody (D01P)

**Catalog Number:** H00147746-D01P

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against a full-length human HIPK4 protein.

**Immunogen:** HIPK4 (NP\_653286.2, 1 a.a. ~ 616 a.a) full-length human protein.

**Sequence:**

MSTIQSETDCYDIEVLGKGTGFEVAKGWRRSTGEMV  
AIKILKNDAYRNRIIKNELKLLHCMRGLDPEEAHVIRFLE  
FFHDALKFYLVFELLEQNLFEFQKENNFAPLPARHIRT  
TLQVLTALARLKLAIHADLKPENIMLVDQTRCPFRVK  
VIDFGSASIFSEVRYVKEPYIQSRFYRAPEILLGLPFCEK  
VDVWSLGCVMAELHLGWPLYPGNNEYDQVRYICETQ  
GLPKPHLLHAACKAHFFKRNPHPDAAANPWQLKSSA  
DYLAETKVRPLERRKYMLKSLDQIETVNGGSVASRLTF  
PDREALAEHADLKSMVELIKRMLTWESHERISPSAALR  
HPFVSMQQLRSAHETTHYYQLSLRSYRLSLQVEGKPP  
TPVVAEDGTPYYCLAEKEAAGMGSVAGSSPFFREE  
KAPGMQRAIDQLDDLSLQEAGHGLWGETCTNAVSDM  
MVPLKAAITGHHVPDSGPEPILAFYSSRLAGRHKARKP  
PAGSKSDSNFSNLIRLSQVSPEDDRPCRGSSWEEGE  
HLGASAEPLAILQRDEDGPNIDNMTMEAERPDPFLFD  
PSSCPGEWLSEPDCTLESVRGPRAQGLPPRRSHQH  
PPRGATSFLQHVTGHH

**Host:** Rabbit

**Reactivity:** Human

**Applications:** WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 147746

**Gene Symbol:** HIPK4

**Gene Alias:** FLJ32818