

Datasheet

ZNF582 purified MaxPab mouse polyclonal antibody (B01P)

Gene Alias: FLJ30927

Catalog Number: H00147948-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF582 protein.

Immunogen: ZNF582 (NP_653291.1, 1 a.a. ~ 517 a.a) full-length human protein.

Sequence:

MSLGSELFDRDVAIVFSQEEWQWLAPAQRDLRYDVML
ETYSNLVSLGLAVSKPDVISFLEQGKEPMMVERVVSG
GLCPVLESRYDTKELFPKQHVYEVESQWEIMESLTS
YGLECSSFQDDWECRNQFDRQQGNPDRHFHQMIIRH
EEMPTFDQHASLTFYQKIHTREKPFQYNGKCRKDFWQK
ELLINHQQGIYTNKPYKCKEKGKAFKYGSRLIQHENIHS
GKKPYECKECKGAFNSGNSFIQHQRVHTGEKPYECK
DCEKAFSRSSQLIEHQRTHTGEKPYQCKECKGAFNRI
SHLKVHYRIHTGEKPYACECKGKTFSHRSQLIQHQTV
HTGRKLYECKECKGAFNQGSTLIRHQRIHTGEKPYEC
KVCGKAFRVSSQLKQHQRHTGEKPYQCKVCGRAFK
RVSHLTVHYRIHTGEKPYECKECKGAFSHCSQLIHHQ
VIHTEKKPYEYKECEKTLSDSTTVQPMHNRETHV
NIINVEKPSISSYPLLIREFMLASNHMNGSNGESPLA

Host: Mouse

Reactivity: Human

Applications: Det Ab, IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 147948

Gene Symbol: ZNF582