

Datasheet

LOC148203 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00148203-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human LOC148203 protein.

Immunogen: LOC148203 (AAH34499.1, 1 a.a. ~ 137 a.a) full-length human protein.

Sequence:

MDDLRYGVYPVKGASGYPGAERNLLEYSYFEKGPLTF
RDVVIEFSQEEWQCLDTAQQDLRKVMLENFRNLVFL
GIDVSKPDLITCLEQKDPWNMKRHSMTATPPESGVL
KFPTIIILLSRCFFQSVNICFIYLEP

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 148203

Gene Symbol: ZNF738

Gene Alias: -