

Datasheet

CREB3L4 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00148327-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human CREB3L4 protein.

Immunogen: CREB3L4 (NP_570968.1, 1 a.a. ~ 395 a.a) full-length human protein.

Sequence:

MDLGIPDLLDAWLEPPEDIFSTGSVLELGLHCPPPEVP
VTRLQEQGLQGWSGGDRGCLQESEPEDFLKLFID
PNEVYCSEASPGSDSGISEDPCHPDSPAPRATSSPM
LYEVVYEAGALERMQGETGPNVGLISIQLDQWSPAFM
VPDSCMVSELPFDAHAHILPRAGTVAPVPCTTLLPCQT
LFLTDEEKRLLGQEGVSLPSHLPLTKAEERVLKKVRRK
IRNKQSAQDSRRRKKEYIDGLESRVAACSAQNQELQK
KVQELERHNISLVAQLRQLQTIAQTSNKAQTSTCVLI
LLFSLALIILPSFSPFQSRPEAGSEDYQPHGVTSRNILT
HKDVTENLETQVVESRLREPPGAKDANGSTRTLLEKM
GGKPRPSGRIRSVLHADEM

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 148327

Gene Symbol: CREB3L4

Gene Alias: AIBZIP, ATCE1, CREB3, CREB4, JAL, hJAL

Gene Summary: cAMP response element-binding (CREB) proteins are a family of mammalian transcription activators. For further background information on CREB proteins, see CREB1 (MIM 123810).[supplied by OMIM]