

Datasheet

CNIH3 monoclonal antibody (M01), clone 1E10

Catalog Number: H00149111-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full length recombinant CNIH3.

Clone Name: 1E10

Immunogen: CNIH3 (AAH22780, 1 a.a. ~ 160 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MAFTFAAFCYMLSLVLCAALIFFAIWHIIAFDELRTDFKS
PIDQCNPVHARERLRNIERICFLLRKLVLPEYSIHSIFCI
MFLCAQEWLTLGLNVPLLFYHFWRYFHCPADSSELAY
DPPVVMNADTLSYCQKEAWCKLAFYLLSFFYYLYCMI
YTLVSS

Host: Mouse

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 149111

Gene Symbol: CNIH3

Gene Alias: FLJ38993