

Datasheet

PDIK1L monoclonal antibody (M01), clone 4B7

Catalog Number: H00149420-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant PDIK1L.

Clone Name: 4B7

Immunogen: PDIK1L (NP_690048, 1 a.a. ~ 99 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MVSSQPKYDLIREVGRGSYGVVYEAVIRKTSARVAVK
KIRCHAPENVELALREFWALSSIKSQHPNVHLEECILQ
KDGMVQKMSHGSNSSLYLQLVET

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 149420

Gene Symbol: PDIK1L

Gene Alias: CLIK1L, RP11-96L14.4