

Datasheet

BNIP1 purified MaxPab mouse polyclonal antibody (B03P)

Catalog Number: H00149428-B03P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human BNIP1 protein.

Immunogen: BNIP1 (AAH74779.3, 1 a.a. ~ 275 a.a) full-length human protein.

Sequence:

MRKRLSAPELRLSLTKGPGNDGASPTQSAPSSPDGSS
DLEIDELETPSDSEQLDSGHEFEWEDELPRAEGLGTS
ETAERLGRGCMWDVTGEDGHHWRVFRMGPREQRV
DMTVIEPYKKVLSHGGYHGDGLNAVILFASCYLPRSSI
PNYTYVMEHLFRYMGVLTLELLVAENYLLVHLSGGTSR
AQVPPLSWIRQCYRTLDRRLRKNLRALVVVHATWYVK
AFLALLRPFISKFTRKIRFLDSLGEAQLISLDQVHIPE
AVRQLDRDLHGSGGT

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 149428

Gene Symbol: BNIP1

Gene Alias: BNIP-S, BNIP1-1, BNIP1-2, BNIP1,
BNIP2, PP753