

Datasheet

LIPI purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00149998-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human LIPI protein.

Immunogen: LIPI (AAI48692.1, 1 a.a. ~ 481 a.a)
full-length human protein.

Sequence:

MLLKCLHNNLCQKYSAHAFQFSPRNVLLVCLRSD
NKRPCLEFSQLSVKDSFRDLFIPRIETILMMYTRNNLNC
AEPLFEQNNSLNVNFNTQKKTVWLIHGYPVGSIPLWL
QNFVRILLNEEDMNVIVVDWSRGATTFIYNRAVKNTRK
VAVSLSVHIKNLLKHGASLDNFHFIGVSLGAHISGFVGK
IFHGQLGRITGLDPAGPRFSRKPPYSRLDYTDKAFVDV
IHSDSNGLGIQEPLGHIDFYPNGGNKQPGCPKSIFSGI
QFIKCNHQRAVHLFMALETNCNFISFPCRSYKDYKTS
LCVDCDCFKKSCPRLGYQAKLFKGVKERMGRPLR
TTVFLDTSGTYPCTYYFVLSIIVPDKTMMDGSSFSKLL
NQLGMIEEPRLYEKNKPFYKLQEVKILAQFYNDFVNISS
IGLTYFQSSNLQCSTCTYKIQRMLKSLTYPERPPLCR
YNIVLKDREEVFLNPNTCTPKNT

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 149998

Gene Symbol: LIPI

Gene Alias: LPDL, PRED5