

Datasheet

SNF1LK purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00150094-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human SNF1LK protein.

Immunogen: SNF1LK (AAH38504.1, 1 a.a. ~ 783 a.a.) full-length human protein.

Sequence:

MVIMSEFSADPAGQGQGGQKPLRVGFYDIERTLGKG
NFAVVKLARHRVTKTQVAIKIIDKTRLDSSNLEKIYREV
QLMKLLNHPHIIKLYQVMETKDMLYIVTEFAKNGEMFD
YLTSGHLSSENEARKKFWQILSAVEYCHDHHIVHRDLK
TENLLLDGNMDIKLADFGFGNFYKSGEPLSTWCGSP
YAAPEVFEGKEYEGPQLDIWSLGVVLYVLVCGSLPFD
GPNLPTLRQRVLEGRFRIPFFMSQDCESLIRRMVLVDP
ARRITIAQIRQHRWMRAEPCLPGPACPAFSAHSYTSNL
GDYDEQALGIMQTLGVDRQRTVESLQNSSYNHFAAIY
YLLERLKEYRNAQCARPGPARQPRPRSSDLSGLEVP
QEGSTDPFRPALLCPQPQTLVQSVLQAEMDCELQSS
LQWPLFFPVDASCSGVFRPRPVSPSSLLDTAISEEAR
QGPGLLEEEDTQESLPSSTGRRHTLAEVSTRLSPLTA
PCIVVSPSTTASPAEGTSSDSCLTFSASKSPAGLSGTP
ATQGLLGACSPVRLASPFLGSQSATPVLQAQGGGLGGA
VLLPVSFQEGRRASDTSLTQGLKAFRQQLRKTRTKG
FLGLNLIKGLARQVCQVPASRASRGGLSPFHAPAQSP
GLHGGAAGSREGWSLLEEVLEQQRLLQLQHHPAAAP
GCSQAPQPAPAPFVIAPCDGPGAAPLPSTLLTSGLPLL
PPPLLQTGASPVASAAQLLDTHLHIGTGPTALPAVPPP
RLARLAPGCEPLGGLLQGDCEMEDLMPCSLGTFVLVQ

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 150094

Gene Symbol: SIK1

Gene Alias: MSK, SIK, SNF1LK